

DENT**OLOGY**

by **ABUTMENT DIRECT**

Same DNA. More Solutions.

ASIGA[®]


AMANN GIRRBACH

keystone[®]
Industries

DFS
Dental Future Systems
DIAMON

BIG CAT
Precision Rotary Instruments

DETAX

HUF SCHMIED
ZERSpanungssysteme

HarvestDental

 **buffalo**

NK-Optik

MESA
PRODUCER OF SPECIAL ALLOYS SINCE 1878

pro**3d**ure
medical

PROFESSIONAL
SNOW ROCK


TALLADIUM, INC.


Whip Mix


MERZ
DENTAL


DEKEMA

4th Edition

Content

3D PRINTING

ASIGA PRINTERS	4
ASIGA PRINT RESIN	8
DETAX PRINT RESIN	6
KEYSTONE PRINT RESIN	8
PRO3DURE PRINT RESIN	9
NK OPTIK LIGHT CURING OVEN	5
ULTRASONIC CLEANER	5
WHIP MIX PRINT RESIN	9
3D PRINTING WASHING FLUIDS	6

CAD/CAM

AMANN GIRRBACH MILLING & SCANNING EQUIPMENT	10
AMANN GIRRBACH SINTERING FURNACES	12
AMANN GIRRBACH MILLING TOOLS	21
HUFSCHMIED MILLING TOOLS	20
DEKEMA FURNACES	13
AMANN GIRRBACH DISCS ZIRCONIA & PMMA	14
HARVEST DENTAL WAX DISCS	18
KEYSTONE KEYMILL DENTURE DISCS	19
MERZ PEEK DISCS	19
TALLADIUM TITANIUM & CHROME-COBALT DISCS	19

LAB SUPPLIES

HARVEST DENTAL LAB MATERIALS	22
KEYSTONE DENTAL LAB MATERIALS	25
DETAX DENTAL LAB MATERIALS	30
DETAX VPS IMPRESSION AND BITE MATERIAL	33
BUFFALO DENTAL LAB MATERIALS	35
SNOW ROCK SCAN POWDER	32
MESA ALLOY	35
TALLADIUM ALLOY	35
BIG CAT FINISHING & POLISHING BURS	41
DFS ROTARY TOOLS	61

LAB EQUIPMENT

AMANN GIRRBACH LAB EQUIPMENT	36
KEYSTONE LAB EQUIPMENT	38
BUFFALO LAB EQUIPMENT	39

Same DNA. More Solutions.

In 2010 Abutment Direct Inc. was created with the mission of supplying high quality and price-effective implant prosthetic components in partnership with Medentika GmbH. As the number of satisfied customers grew across Canada we progressively expanded our product solutions outside of the core implant component market.

The time is right to introduce even more solutions by expanding our product offerings and doing so through our new division, Dentology. What remains the same is our goal of working with premier global manufacturing partners, and our commitment to fast, friendly and knowledgeable customer service.

Mike Parsons
President
Abutment Direct Inc.

Follow us:



@dentologybyabutmentdirect

Terms and Conditions of Sales

Pricing:

- Pricing is subject to change without notice. Please contact us to confirm current pricing.
- Pricing excludes applicable taxes.

Shipping:

- Shipping is free on orders over \$500 before taxes.
- Please contact us for current standard shipping charges to your location.
- Additional shipping charges may apply for bulky or heavy items.

Returns and credits:

- Items that are unopened, unused and undamaged and still in the original packaging may be returned for a full credit to the customer's account within six weeks of purchase date, less a 10% restocking fee.
- Packaging must be in its original condition, with no markings added by the customer. All credits will be subject to inspection of the goods upon their return.
- Goods are returned at the customer's expense using the secure method of their choice. If goods have been sent to the customer in error Dentology will arrange for pick up and return of the product at Dentology's expense.

ASIGA®

Printers

The #1 Selling 3D Printers in Canada

MAX UV Printer



The ultimate 3D DLP printer for dental professionals needing to produce a maximum of approximately 20 pieces per day **Call for pricing.**

Product Specification	
Build Volume X, Y, Z	119 × 67 × 75mm (4.69 × 2.64 × 2.95 inches)
Pixel Resolution	62µm
Technology	DLP
LED Wavelength	385nm (high power UV LED)
Material Compatibility	Asiga materials, 100% Open Material System - use any suitable third party material
Software	Asiga Composer software. Lifetime updates included
File Inputs	STL, SLC, PLY, STM (Asiga Stomp file format)
Network Compatibility	Wifi, Wireless Direct, Ethernet
Power Requirements	100 to 240 VAC, 50/60 Hz, ≤2 Amps
System Sizing	260 × 385 × 370mm / 19.3 kg (10.2 × 15.2 × 14.5 inches / 42.6 lb)
Packed Sizing	410 × 500 × 480mm / 21.6 kg (16.1 × 19.7 × 18.9 inches / 47.6 lb)
Warranty	12 months manufacturer's warranty
Technical Support	Unlimited lifetime technical support included
Bundle Includes	3D printer, Composer software, 1 kg Asiga material, 1L build tray, calibration toolkit

PRO 4K

Larger scale 3D DLP printer for dental professionals needing to produce greater than 20 pieces per day. **Call for pricing.**



PRO 4K80 Product Specification	
Build Volume X, Y, Z	217 × 122 × 200mm (8.54 × 4.8 × 7.87 inches)
Pixel size – 4K mode	56µm
Pixel size – Native Mode	80µm
Technology	DLP with pixel-shifting technology
LED Wavelength	385nm (high power UV LED) or 405nm
Material Compatibility	Asiga materials, 100% Open Material System — use any suitable third party material
Software	Asiga Composer software. Lifetime updates included
File Inputs	STL, SLC, PLY, STM (Asiga Stomp file format)
Network Compatibility	Wifi, Wireless Direct, Ethernet
Power Requirements	100-240VAC, 50/60Hz, 500 Watts (100V - 5Amp Max. 240V - 2.1Amp)
System Sizing	465 × 420 × 1370mm / 140 kg 18.3 × 16.5 × 53.9 inches / 309 lb
Packed Sizing	900 × 700 × 1540mm / 205 kg 35.4 × 27.6 × 60.6 inches / 452 lb
Warranty	12 months manufacturer's warranty
Technical Support	Unlimited lifetime technical support included
Bundle Includes	3D printer, Composer software, 1 kg Asiga material, 1L build tray, calibration toolkit

Printer Accessories

Item	Description	Printer	Price
Resin Tray	1Litre	MAX UV	\$125
Resin Tray	2Litre	MAX UV	\$180
Resin Tray	5Litre	MAX UV	\$300
Resin Tray	10Litre	MAX UV	\$400
Resin Tray Storage		MAX UV	\$19
Resin Tray	2Litre	PRO 4K	\$125
Resin Tray	5Litre	PRO 4K	\$180
Resin Tray	10Litre	PRO 4K	\$300
Resin Tray Storage		PRO 4K	\$49
Build Plate		MAX UV	\$699
Build Plate		PRO 4K	\$ Call



Otoflash G171

Flash-curing unit for printed parts

CALL FOR PRICING

- Two flash bulbs produce ten flashes per second
- Light radiation at a wavelength of 280-580 nm
- Average flash power 200W
- Suitable for curing all light-cured materials
- Option for nitrogen gas injection to eliminate oxygen inhibition layer
- Includes gas regulator
- Dimensions: 310 x 310 x 140 mm
- Weight: 6 kg.
- Replacement tray
- Filter
- Replacement bulb module



Ultrasonics for Printed Part Cleaning

Dentology JPS-30A 6 Litre Ultrasonic Cleaner

CALL FOR PRICING

- Stainless Steel construction
- Heating Power: 200W
- Ultrasonic Power: 180W
- Frequency: 40kHz
- Power Supply: AC110V/60Hz
- Time Setting: 1 - 30 minutes (LED Digital Display)
- Temperature Setting: 0-80 degrees Celsius
- Transducer Quantity: 3 set
- Tank Size: 11 3/4" x 6" x 6" (L x W x H)
- Ultrasonic Cleaner Size: 14 1/2" x 6 3/4" x 10 3/4" (L x W x H)





Optimize 99% IPA
Applications:
 • Isopropyl alcohol for printed part cleaning


Quantity: 943ml
\$19



Innovaprint Wash 3D
Applications:
 • Green alternative cleaning solution for printed parts
 • No IPA
 • No sell or fumes
 • Resusable and biodegradable


Quantity: 1L, 5L
 1L **\$39**
 5L **\$139**

DETAX 3D Print Resins




FREEPRINT[®] denture
 Resin for the 3D printing of removable denture bases. Colour: pink-transparent

Product Number and Name	Unit	Price
02060 FREEPRINT denture pink-transparent	500g	\$259
02040 FREEPRINT denture pink-transparent	1kg	\$489
04092 FREEPRINT denture pink	1kg	\$489




FREEPRINT[®] temp
 Resin for the 3D printing of temporary crowns & bridges. Shades: A1, A2, A3

Product Number and Name	Unit	Price
04058 FREEPRINT temp 385 A1	500g	\$319
04059 FREEPRINT temp 385 A2	500g	\$319
04060 FREEPRINT temp 385 A3	500g	\$319
04062 FREEPRINT temp 385 A1	1kg	\$549
04063 FREEPRINT temp 385 A2	1kg	\$549
04064 FREEPRINT temp 385 A3	1kg	\$549




FREEPRINT[®] Splint 2.0
 Resin for the 3D printing of splints, fixation and transfer keys. Colour: clear-transparent

Product Number and Name	Unit	Price
02076 Freeprint Splint 2.0	1kg	\$369




FREEPRINT[®] ortho
 Resin for the 3D printing of base parts for orthodontic apparatuses, drilling and X-ray templates, occlusal and fixation splints. Heat sterilizable. Colour: clear-transparent

Product Number and Name	Unit	Price
03989 FREEPRINT ortho 385	1kg	\$379



FREEPRINT[®] tray
 Resin for the 3D printing of individual impression and functional trays and base resin plates. Colour: green

Product Number and Name	Unit	Price
04086 FREEPRINT tray 385	1kg	\$349



FREEPRINT[®] cast 2.0
 Resin for the 3D printing of high-precision casting objects. Colour: red

Product Number and Name	Unit	Price
02548 FREEPRINT cast 2.0 385	500g	\$149
02632 FREEPRINT cast 2.0 385	1kg	\$269


DETAX 3D Print Resins



FREEPRINT® gingiva

Resin for the 3D printing of flexible gingival masks for dental models.
Colour: gingiva

Product Number and Name	Unit	Price
02820 FREEPRINT gingiva 385	500g	\$149
02843 FREEPRINT gingiva 385	1kg	\$259



FREEPRINT® model 2.0

Resin for the 3D printing of dental models.
Colours: sand, light-grey, grey, caramel, white

Product Number and Name	Unit	Price
02128 FREEPRINT model 2.0 385 sand	1kg	\$229
02099 FREEPRINT model 2.0 385 light-grey	1kg	\$229
02177 FREEPRINT model 2.0 385 grey	1kg	\$229
02850 FREEPRINT model 2.0 385 caramel	1kg	\$229
02148 FREEPRINT model 2.0 385 white	1kg	\$229



FREEPRINT® model T

Temperature-resistant (up to 200C) resin for the 3D printing of dental models for the thermoforming technique.
Colour: light blue

Product Number and Name	Unit	Price
02332 FREEPRINT model T 385	1kg	\$229

ASIGA®

Printer Resins



DentaMODEL

Colour: Almond

Applications:

- Dental models
- Models with removable dies
- Models with forming aligners

Printer compatibility: 385nm and 405nm

Quantity: 1kg

\$245



DentaTRY

Colour: Available in A1, A2, A3, B1, B2, B3

Applications:

- Denture Try-ins
- Custom Impression trays

Printer compatibility: 385nm and 405nm

Quantity: 1kg

\$245



DentaGUM

Colour: Pink

Applications:

- Gingiva mask

Printer compatibility: 385nm and 405nm

Quantity: 1kg

\$245



DentaCAST

Colour: Green

Applications:

- Partial frameworks
- C&B pressing / casting

Printer compatibility: 385nm and 405nm

Quantity: 1kg

\$245

More options available on the next page.

ASIGA®

Printer Resins



DentaFORM

Colour: Grey
Applications:

- Anti-stick formulation ideal for vacuum thermoforming
- Lower cost model production when less precision required

Printer compatibility:
385nm and 405nm
Quantity: 1kg
\$165



DentaSTUDY

Colour: White
Applications:

- Diagnostic and orthodontic study models

Printer compatibility:
385nm and 405nm
Quantity: 1kg
\$245



DentaTRAY

Colour: Turquoise
Applications:

- Custom impression trays

Printer compatibility:
385nm and 405nm
Quantity: 1kg
\$245

keystone® keyprint®

Industries



Product Name	Unit	Price
KeyCast castable resin	1kg	\$329

Best Seller



Product Name	Unit	Price
KeySplint soft splint resin	1kg	\$529



Product Name	Unit	Price
KeyMask soft tissue resin	500g	\$239



Product Name	Unit	Price
KeyTray custom tray resin	1kg	\$289



Product Name	Unit	Price
KeyOrthoModel model resin, Grey	1kg	\$245



Product Name	Unit	Price
KeyGuide surgical guide resin	1kg	\$319



Product Name	Unit	Price
Keymodel Ultra, ivory and light grey	1kg	\$119
	4kg	\$439



Product Name	Unit	Price
KeySplint Hard	1kg	\$529



Whip Mix Resins



Verimodel OS
Colour: Ivory, Golden Brown, White
Printer compatibility: 385nm and 405nm
Quantity: 1kg
\$259



VeriSplint OS
Colour: Clear
Printer compatibility: 385nm
Quantity: 1kg
\$525



Veritray OS
Colour: Neon Green
Printer compatibility: 385nm
Quantity: 1kg
\$359

Scanning

ceramill® map 600


Ceramill map 600 lab scanner - **Call for quote**

The fully automatic high-performance scanner for open articulator scanning

- Integrated, universal fixator for all common types of articulators ensures maximum flexibility
- Intelligent scan height control for optimum model alignment in the scan field offers the highest possible process reliability
- HD scan via 3D sensor with blue light technology and variable resolution guarantees optimal and reproducible results
- DNA Speed Scanning enables a full jaw scan with unrestrictedly usable results in only 18 seconds

Technical Data:

- Axes: 3
- Dimensions: 415 x 424 x 469 mm
- Weight: 25 kg
- Accuracy: 4 µm
- Accuracy (µm) implant bars: 10
- Full arch scan (sec.): 18
- Single die with Ceramill Multi Cap (sec.): 4.5



Milling

ceramill® motion 2

Ceramill motion 2 milling machine - **Call for quote**

AVAILABLE IN DRY AND HYBRID (WET/DRY) VERSIONS

- High-performance milling, grinding, carving thanks to DNA milling/grinding strategies
- Maximum range of materials and indications thanks to 5-axis wet and dry processing in one unit
- Innovative processing procedures ensure maximum value creation (e.g. rotational milling of titanium, processing of full dentures)
- Intelligent machine design guarantees optimum protection of all electronic components during wet operation

Technical Data:

- **Axes:** 5
- **Processing mode:** Dry only or Hybrid (wet/dry)
- **Tool management:** 6-compartment tool changer, management of replacement tools
- **Dimensions:** 595 x 530 x 780 (mm)
- **Weight:** 78kg
- **Power Output:** 750 W
- **Motor Speed:** 100,000 rpm
- **Compressed Air:** 50 L/min. at 6 bar



ceramill® motion 3

Ceramill Motion 3 - **Call for quote**

The most intelligent digital 5-axis single disc hybrid mill in the World

Main features:

- Direct user interface
- Autocalibration
- HD onboard camera
- C-Clamp option

Technical Data:

- **Axes:** 5
- **Processing mode:** Hybrid (wet/dry)
- **Tool management:** 6-compartment tool changer, RFID management of replacement tools
- **Dimensions:** 595 x 530 x 780 (mm)
- **Weight:** 78kg
- **Power Output:** 750 W
- **Motor Speed:** 100,000 rpm
- **Compressed Air:** 50 L/min. at 6 bar





ceramill® matik

Ceramill matik milling machine - Call for quote

- Maximum time savings in resource management through automatic management of tools and materials
- Highest possible convenience due to autonomous wet and dry change and automatic cleaning
- Full range of applications due to maximum indication and material portfolio
- Future-safe due to constant software development and expandable holder portfolios



Technical Data:

- **Axes:** (processing + management): 5 + 4
- **Processing mode:** autonomous dry/wet
- **Spindle torque:** 11.2 Ncm
- **Spindle speed:** 100.000min⁻¹
- **Tool management:** replaceable 27x tool stock trays with RFID management, management of replacement tool
- **Material management:** automatic material management with RFID-supported material holders, 36-fold blank tank
- **Dimensions:** 780 x 1265 x 2000mm
- **Weight of system (empty/max):** 520/625k
- **Compressed air peak:** 200l/min at 6bar in the cleaning mode

Sintering



ceramil® therm 3

Ceramil Therm 3 sintering furnace - Call for quote

Features:

- 250 sintering program positions of which 3 are validated AG programs
- 3 stackable sintering bowls for maximum utilisation of the furnace
- Sintering at the press of a button – very easy operation with touch-screen technology
- USB interface for program updates and log files for monitoring firing
- HD onboard camera
- C-Clamp option

Technical Data:

- **Dimensions:** 500 x 461 x 480 (mm)
- **Weight:** 30kg
- **Electrical connected value:** V/Hz 220- 240/50-60
- **Power Output:** 3.5 kW



ceramil®

Ceramil Argotherm 2 sintering furnace for Sintron - Call for quote

Features:

- Sintering program developed specifically for Ceramil Sintron
- A validated system guarantees consistent high quality of the results
- Active cooling after the sintering process
- Minimal consumption of protective gas
- Sintering at the push of a button – extremely easy operation with touch screen technology
- USB interface for program updates and log files for monitoring firing

Technical Data:

- **Dimensions:** 468 x 461 x 480 (mm)
- **Weight:** 30kg
- **Electrical connected value:** V/Hz 220- 240/50-60
- **Power Output:** 3.5 kW



ceramil®

Ceramil Therm DRS high-speed sintering furnace - Call for quote

Features:

- High-Speed sintering cycle of restorations with Zolid DRS in approx. 20 minutes without any relevant loss of esthetics and strength
- Innovative operating concept with individual sintering programs provides high comfort and safe workflows
- Ideally suited for the rapid provision of single-tooth restorations through perfect workflow integration

Technical Data:

- **Dimensions:** 390 x 300 x 484 (mm)
- **Weight:** 50kg
- **Electrical connected value:** 230V / 50-60Hz
- **Power Output:** 3.5 kW
- **Maximum temperature:** 1,600 °C



In addition to the usual benefits of the standard line of Dekema furnaces the Austromat i-series furnaces feature:

- 10.4" multi-touch display for ease of use
- Internet/LAN connection allowing access to continuous program updates and custom programs, and easy remote access




Austromat 664i sintering furnace - Call for quote

- Maximum temperature 1560C
- Speed mode
- Space for up to 6 unit bridge or up to 20 units in one level – optional accessory to sinter on two levels
- 65mm diameter chamber
- 230V
- Ideal for lower sintering volume requirements



Austromat 674i sintering furnace - Call for quote

- Maximum temperature 1560C
- Speed mode
- Space for full arch bridge or up to 40 units in one level – optional accessory to sinter on two levels
- 100mm diameter chamber
- 230V
- Ideal when manufacturing full arch bridges or for higher volume sintering requirements



Austromat 654i combination firing and pressing furnace - Call for quote

- Maximum temperature 1200C
- Press, fire and glaze with or without vacuum
- Crystallize and infiltrate with or without vacuum
- Press multiple plungers at once with optional trixpress feature
- Certified to press PEEK materials PEKTON and BioHPP
- 230V
- For the lab looking for a versatile, dual-function furnace (press or fire)



FERRO Calibration Ring

Part # PTCR-HTH

Use to calibrate sintering furnace firing temperature

- 1450-1750C

\$12 each

\$99 for ten-pack



Milling Materials



SUPER HIGH TRANSLUCENT Zirconia (SHT)



Zolid fx white

Zolid FX White. Super high translucent zirconia blanks used for monolithic or anatomically reduced anterior restorations and up to 3-unit bridges in the molar region.

- Available in 98.5 mm diameter discs
- Flexural strength 700 Mpa

Part number	Thickness	Price
760130	12mm	\$169
760131	14mm	\$179
760132	16mm	\$189
760133	18mm	\$199
760134	20mm	\$209



Zolid fx preshades

Zolid FX Preshades are pre-stained, super high translucent zirconia blanks for monolithic or anatomically reduced anterior restorations and up to 3-unit bridges in the molar region.

- Available in 98.5 mm diameter discs
- Shades: bleach, light A, medium A, light B, light C and light D
- Flexural strength 700 Mpa



Part number	Thickness	Shade	Price
761551	14mm	Bleach	\$179
761552	16mm	Bleach	\$189
761554	20mm	Bleach	\$199
761556	14mm	A Light	\$179
761557	16mm	A Light	\$189
761559	20mm	A Light	\$199
761561	14mm	A Medium	\$179
761562	16mm	A Medium	\$189
761564	20mm	A Medium	\$199

Part number	Thickness	Shade	Price
761566	14mm	B Light	\$179
761567	16mm	B Light	\$189
761569	20mm	B Light	\$199
761576	14mm	C Light	\$179
761577	16mm	C Light	\$189
761579	20mm	C Light	\$199
761586	14mm	D Light	\$179
761587	16mm	D Light	\$189
761589	20mm	D Light	\$199

Zolid fx multilayer

Zolid FX Multilayer is a polychromatic, super high translucent zirconia with integrated shade and translucency gradient. Smooth shade transitions simulate tooth enamel, dentine and cervical shades in a naturally fluent progression without disruptive shade edges.

- Available in 98.5 mm diameter discs
- Shades O/A1, A2/A3, A3.5/A4, O/B1, B2/B3, B3/B4, C1/C2, C3/C4, D2/D3, D3/D4
- Flexural strength 700 Mpa



Part number	Thickness	Shade	Price
761776	14mm	O/A1	\$189
761777	16mm	O/A1	\$209
761778	20mm	O/A1	\$239
761780	14mm	A2/A3	\$189
761781	16mm	A2/A3	\$209
761782	20mm	A2/A3	\$239
761784	14mm	A3.5/A4	\$189
761785	16mm	A3.5/A4	\$209
761786	20mm	A3.5/A4	\$239
761788	14mm	O/B1	\$189

Part number	Thickness	Shade	Price
761789	16mm	O/B1	\$209
761790	20mm	O/B1	\$239
761792	14mm	B2/B3	\$189
761793	16mm	B2/B3	\$209
761794	20mm	B2/B3	\$239
761796	14mm	B3/B4	\$189
761797	16mm	B3/B4	\$209
761798	20mm	B3/B4	\$239
761800	14mm	C1/C2	\$189
761801	16mm	C1/C2	\$209

Part number	Thickness	Shade	Price
761802	20mm	C1/C2	\$239
761804	14mm	C3/C4	\$189
761805	16mm	C3/C4	\$209
761806	20mm	C3/C4	\$239
761808	14mm	D2/D3	\$189
761809	16mm	D2/D3	\$209
761810	20mm	D2/D3	\$239
761812	14mm	D3/D4	\$189
761813	16mm	D3/D4	\$209
761814	20mm	D3/D4	\$239



HIGH TRANSLUCENT Zirconia (HT)

Zolid Gen-X Multi-Layer zirconia

NEW

The new Zolid Gen-X is a multilayer, highly translucent zirconia and is a true all-rounder in terms of "universal use". Due to its high strength and natural translucency, this material can be used for virtually every indication.

- Available in 98.5mm discs
- Flexural strength 1,000 MPa
- Available in shades BL3, BL1, A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4



Part number	Thickness	Shade	Price
766809	14mm	BL3	\$199
766810	16mm	BL3	\$219
766812	20mm	BL3	\$249
766814	25mm	BL3	\$289
766791	14mm	BL1	\$199
766792	16mm	BL1	\$219
766794	20mm	BL1	\$249
766796	25mm	BL1	\$289
766818	14mm	A1	\$199
766819	16mm	A1	\$219
766821	20mm	A1	\$249
766823	25mm	A1	\$289
766827	14mm	A2	\$199
766828	16mm	A2	\$219
766830	20mm	A2	\$249
766832	25mm	A2	\$289
766836	14mm	A3	\$199
766837	16mm	A3	\$219
766839	20mm	A3	\$249
766841	25mm	A3	\$289
766845	14mm	A3.5	\$199
766846	16mm	A3.5	\$219
766848	20mm	A3.5	\$249
766859	25mm	A3.5	\$289

Part number	Thickness	Shade	Price
766854	14mm	A4	\$199
766855	16mm	A4	\$219
766857	20mm	A4	\$249
766859	25mm	A4	\$289
766863	14mm	B1	\$199
766864	16mm	B1	\$219
766866	20mm	B1	\$249
766868	25mm	B1	\$289
766872	14mm	B2	\$199
766873	16mm	B2	\$219
766875	20mm	B2	\$249
766877	25mm	B2	\$289
766881	14mm	B3	\$199
766882	16mm	B3	\$219
766884	20mm	B3	\$249
766886	25mm	B3	\$289
766890	14mm	B4	\$199
766891	16mm	B4	\$219
766893	20mm	B4	\$249
766895	25mm	B4	\$289
766899	14mm	C1	\$199
766900	16mm	C1	\$219
766902	20mm	C1	\$249
766904	25mm	C1	\$289

Part number	Thickness	Shade	Price
766908	14mm	C2	\$199
766909	16mm	C2	\$219
766911	20mm	C2	\$249
766913	25mm	C2	\$289
766917	14mm	C3	\$199
766918	16mm	C3	\$219
766920	20mm	C3	\$249
766922	25mm	C3	\$289
766926	14mm	C4	\$199
766927	16mm	C4	\$219
766929	20mm	C4	\$249
766931	25mm	C4	\$289
766935	14mm	D2	\$199
766936	16mm	D2	\$219
766938	20mm	D2	\$249
766940	25mm	D2	\$289
766944	14mm	D3	\$199
766945	16mm	D3	\$219
766947	20mm	D3	\$249
766949	25mm	D3	\$289
766953	14mm	D4	\$199
766954	16mm	D4	\$219
766956	20mm	D4	\$249
766958	25mm	D4	\$289

Milling Materials



Zolid HT+ white

Zolid HT+ is a highly translucent zirconium oxide which combines high mechanical characteristics with outstanding esthetics.

- Available in 98.5 mm diameter discs
- Flexural strength 1100 Mpa

Part number	Thickness	Price
766020	10mm	\$139
766021	12mm	\$149
766022	14mm	\$169
766023	16mm	\$179
766025	20mm	\$189
766027	25mm	\$209



HIGH TRANSLUCENT Zirconia (HT)

Zolid HT+ preshades

- Available in 16 Vita A-D shades and 98.5mm discs
- Flexural strength 1000 Mpa



Part number	Thickness	Shade	Price
766152	14mm	Bleach	\$169
766153	16mm	Bleach	\$179
766155	20mm	Bleach	\$199
766159	14mm	A1	\$169
766160	16mm	A1	\$179
766162	20mm	A1	\$199
766166	14mm	A2	\$169
766167	16mm	A2	\$179
766169	20mm	A2	\$199
766173	14mm	A3	\$169
766174	16mm	A3	\$179
766176	20mm	A3	\$199
766180	14mm	A3.5	\$169
766181	16mm	A3.5	\$179
766183	20mm	A3.5	\$199
766187	14mm	A4	\$159
766188	16mm	A4	\$169

Part number	Thickness	Shade	Price
766190	20mm	A4	\$199
766194	14mm	B1	\$169
766195	16mm	B1	\$179
766197	20mm	B1	\$199
766201	14mm	B2	\$169
766202	16mm	B2	\$179
766204	20mm	B2	\$199
766208	14mm	B3	\$169
766209	16mm	B3	\$179
766211	20mm	B3	\$199
766215	14mm	B4	\$169
766216	16mm	B4	\$179
766218	20mm	B4	\$199
766222	14mm	C1	\$169
766223	16mm	C1	\$179
766225	20mm	C1	\$199
766229	14mm	C2	\$169

Part number	Thickness	Shade	Price
766230	16mm	C2	\$179
766232	20mm	C2	\$199
766236	14mm	C3	\$169
766237	16mm	C3	\$179
766239	20mm	C3	\$199
766243	14mm	C4	\$169
766244	16mm	C4	\$179
766246	20mm	C4	\$199
766250	14mm	D2	\$169
766251	16mm	D2	\$179
766253	20mm	D2	\$199
766257	14mm	D3	\$169
766258	16mm	D3	\$179
766260	20mm	D3	\$199
766264	14mm	D4	\$169
766265	16mm	D4	\$179
766267	20mm	D4	\$199



LOW TRANSLUCENT Zirconia

Zolid ZI

Lowest translucency zirconia for use in fabricating abutments, to block out dark stump shades and to mask metals

- Available in 98.5 mm diameter discs
- Flexural strength 1200 Mpa

Part number	Thickness	Price
760221	14mm	\$154
760223	16mm	\$164
760225	18mm	\$174
760227	20mm	\$184
760229	25mm	\$199



TEMP AND SPLINT

Ceramill A-Temp PMMA discs



- Monochromatic
- 98.5 mm
- Shades A1, A2, A3, A3.5, B2, C2
- Flexural strength 135 Mpa

Part number	Thickness	Shade	Price	Part number	Thickness	Shade	Price
761350	14mm	A1	\$85	761359	14mm	A3.5	\$85
761351	16mm	A1	\$95	761360	16mm	A3.5	\$95
761352	20mm	A1	\$105	761361	20mm	A3.5	\$105
761353	14mm	A2	\$85	761362	14mm	B2	\$85
761354	16mm	A2	\$95	761363	16mm	B2	\$95
761355	20mm	A2	\$105	761364	20mm	B2	\$105
761356	14mm	A3	\$85	761365	14mm	C2	\$85
761357	16mm	A3	\$95	761366	16mm	C2	\$95
761358	20mm	A3	\$105	761367	20mm	C2	\$105



Ceramill A-Temp Multilayer PMMA discs



- Multilayer
- 98.5 mm
- Shades O/A1, A2/A3, B2/B3, C1/C2
- Flexural strength 135 Mpa

Part number	Thickness	Shade	Price	Part number	Thickness	Shade	Price
761055	14mm	O/A1	\$119	761061	14mm	B2/B3	\$119
761056	16mm	O/A1	\$129	761062	16mm	B2/B3	\$129
761057	20mm	O/A1	\$149	761063	20mm	B2/B3	\$149
761058	14mm	A2/A3	\$119	761064	14mm	C1/C2	\$119
761059	16mm	A2/A3	\$129	761065	16mm	C1/C2	\$129
761060	20mm	A2/A3	\$149	761066	20mm	C1/C2	\$149



Ceramill A- Splint Transparent PMMA discs



- Transparent PMMA for splints and try-ins
- 98.5mm
- Flexural strength >100 Mpa

Part number	Thickness	Price
761084	14mm	\$59
761085	16mm	\$69
761086	20mm	\$89



Milling Materials

Ceramill Sintron Co-Cr Discs - easy as wax

Ceramill Sintron revolutionizes CNC-based dry processing of non-precious metal restorations using desktop milling machines in-house in the laboratory. The wax-like texture of Ceramill Sintron allows the material to be dry milled effortlessly. During the subsequent sintering process under shielding gas atmosphere the frameworks achieve their final material properties – a CoCr unit with very homogeneous material structure. Any conventional CoCr framework porcelain can then be used for veneering. Ceramill Sintron was developed in collaboration with the Fraunhofer IFAM Dresden, Germany.

- Available in 98.5mm discs
- 66% Co, 28% Cr, 5% Mo
- Tensile strength 900 MPa



Part Number	Thickenss	Price
761132	10mm	\$369
761130	12mm	\$399
761128	14mm	\$429
761126	16mm	\$459
761124	18mm	\$499
761122	20mm	\$599
761120	25mm	\$699

Wax Press 98 pressing wax

Millable wax for pressing and casting. Includes option for addition wax to repair imperfections.



Size	Beige	Grey	White	Price
98.5 x 12mm	10002-12	10004-12	10007-12	\$39
98.5 x 14mm	10002-14	10004-14	10007-14	\$39
98.5 x 16mm	10002-16	10004-16	10007-16	\$39
98.5 x 18mm	10002-18	10004-18	10007-18	\$39
98.5 x 20mm	10002-20	10004-20	10007-20	\$39
3 x 20g Addition Wax tins	20850	20852	20856	\$45
1 x 20g Addition Wax tin	20820	20822	20826	\$29

Wax press Esthetic 98 esthetic pressing wax

For beautiful and realistic looking aesthetic wax ups.

- Dentin shaded
- Can be cutback and refined with addition wax
- Pressable



Size	BL-2	Clear	Price
98.5 x 12mm	10010-12	N/A	\$39
98.5 x 14mm	10010-14	N/A	\$39
98.5 x 16mm	10010-16	N/A	\$39
98.5 x 18mm	10010-18	N/A	\$39
98.5 x 20mm	10010-20	N/A	\$39
1 x 20g Addition Wax tin	-----	20826	\$22
3 x 20g Addition Wax tins	20859	20856	\$45

Purity Wax 98 carvable pressing wax

Unique millable and carvable pressing wax that allows the digital technician to easily refine the milled output using wax carving instruments.

- Precise milling.
- Easy carving.
- Perfect adaptability.



Size	Beige	Price
98.5 x 12mm	15001-12	\$49
98.5 x 14mm	15001-14	\$49
98.5 x 16mm	15001-16	\$49
98.5 x 18mm	15001-18	\$49
98.5 x 20mm	15001-20	\$49
98.5 x 25mm	15001-25	\$49

Milling Materials

keymill™
by keystone industries

Denture base disc

KeyMill high-impact denture milling disc, featuring Diamond D.



Size	Product	Price
25mm	Original or Light Reddish pink	\$79
30mm	Original or Light Reddish pink	\$92



PEEK Biosolution discs

- Flexural strength – 170 MPa (ISO-20795-1)
- Elastic Modulus – 5,100 MPa (white and A2/B2), 4,500 MPa (Gum)
- PEEK and titanium dioxide



*S = Stock
SO = Special Order

Product Number	Thickness	Shade	Type*	Price
1018960	10mm	white	SO	\$389
1018961	12mm	white	SO	\$409
1018962	16mm	white	S	\$449
1018963	20mm	white	S	\$499
1018964	24mm	white	S	\$539
1018984	26mm	white	SO	\$579
1018985	30mm	white	SO	\$639
1018965	10mm	A2/B2	SO	\$389
1018966	12mm	A2/B2	SO	\$409
1018967	16mm	A2/B2	S	\$449
1018968	20mm	A2/B2	S	\$499
1018969	24mm	A2/B2	S	\$539
1018990	26mm	A2/B2	SO	\$579
1018991	30mm	A2/B2	SO	\$639
1019230	16mm	Gum	S	\$449
1019231	20mm	Gum	S	\$499
1019232	24mm	Gum	S	\$539
1019234	30mm	Gum	SO	\$639

Luminesse Metal Milling Discs




	Thickness	Chrome Cobalt	Ti-Cam Grade 5
98.5mm discs	10mm	\$279	\$119
	12mm	\$299	\$149
	14mm	\$309	\$179
	16mm	\$339	\$219
	18mm	\$359	\$259
	20mm	\$369	n/a

Milling Machine Tools

HUFSCHMIED

ZERSpanungSSYSTEME



- Precision-manufactured in Germany
- Specialty proprietary coatings for extended tool life
- DPI coating for zirconia is 100% diamond for superior surface finish
- All cutters **\$129 each**

	Product Number	Dimensions				
		D1 mm	D2 mm	L2 +1 mm	L1 + 2mm	L3 + 1mm
AMANN GIRRBACH						
CeraMill Mikro and Motion 2						
PMMA/WAX/Peek (Ball end mill 1F)						
	110FPRDTA006-6-AG	0.6	3	0.8	47	6
	110FPRDTAL01-16-AG	1	3	2	47	16
	110FPRDTAL025-16-AG	2.5	3	5	47	6
Zirconia (Ball end mill 2F) Diamond coated						
	HM915ECOK003-3-AG-DIP	0.3	3	0.4	47	6
	HM915ECOK006-6-AG-DIP	0.6	3	0.6	47	10
	HM915VL010-AG-DIP	1	3	1	47	16
	HM915L025-AG-DIP	2.5	3	3	47	16
IMES-ICORE						
Coritec 350i, 550i, 650i, 850i						
PMMA/WAX/Peek (Ball end mill 1F)						
	110-PRDTA006060-6-ARM	0.6	6	0.8	53	6
	110-PRDTA01060-13-ARM	1	6	2	53	13
	110-PRDTA015060-18-ARM	1.5	6	3	53	18
	110-PRDTA0250060-18-ARM	2.5	6	5	53	18
Cobalt-Chrome (Ball end mill 2F) TL 12 coated						
	HC632DTA010060-10-ARM	1	6	2	53	10
	HC632DTA015060-10-ARM	1.5	6	3	53	10
	HC632DTA020060-16-ARM	2	6	4	53	16
	HC632DTA025060-16-ARM	2.5	6	5	53	16
Zirconia (Ball end mill 2F) Diamond coated						
	HM915DTMFEL005-ARM-DIP	0.6	6	0.6	53	4.5
	HM915DTMFEL010-ARM-DIP	1	6	1	53	8
	HM915DTMFEL015-ARM-DIP	1.5	6	1.5	53	9.5
	HM915DTMFXL025-ARM-DIP	2.5	6	3	53	20
ROLAND						
DWX-50, DWX-51D, DWX-52D, DWX-52DCi						
PMMA/WAX/Peek (Ball end mill 1F)						
	110FPRDTA006040-6	0.6	4	0.8	50	6
	110FPRDTA010040-15	1	4	2	50	15
	110FPRDTA020040-20	2	4	4	50	20
Zirconia (Ball end mill 2F) Diamond coated						
	HM915DTMF003-DIP	0.3	4	0.3	50	3
	HM915DTMFL006-DIP	0.6	4	0.6	50	12
	HM915DTMFVL010-DIP	1	4	1	50	15
	HM915DT020040-20-DIP	2	4	2	50	20
ZIRKONZAHN						
M1 soft						
PMMA/WAX/Peek (Ball end mill 2F)						
	110FRDTAXL003-ZZ	0.3	3	3	57	6
	110FRDTAXL005-ZZ	0.5	3	3	57	7.4
	110FRDTAXL010-ZZ	1	3	4.5	57	7.4
	110FRDTAXL020-ZZ	1.5	3	7.5	57	17
Zirconia (Ball end mill 2F) Diamond coated						
	HM915XL003030-ZZ-DIP	0.3	3	0.3	57	6
	HM915XL005030-ZZ-DIP	0.5	3	0.8	57	7.5
	HM915XL010030-ZZ-DIP	1	3	6	57	12.3
	HM915XL020030-ZZ-DIP	1.5	3	10	57	18

Tool Matrix



	Art. Nr.	Description	Ring color	Ceramill DNA Milling Units						Material															
				Mikro 4X	Mikro 5X	Mikro IC	Motion 2	Matik	Zolid Zirconia	Ceramill Sintron	Ceramill M-Plast	Ceramill PMMA	Ceramill Peek	Ceramill Wax	Ceramill D-Set	Ceramill Splintec	Ceramill D-Wax	Ceramill TI-Forms	Ceramill COCR	Glaskeramik	Lithium Disilikate	Hydrieramik	VITA VIONIC/BDS	Solvay Ultraire AKP	VITA VIONIC Wax White
	760660	Roto RFID 2,5 ZI	grey	x	x		x	x	x																
	760661	Roto RFID 1,0 ZI	grey	x	x		x	x	x																
	760662	Roto RFID 0,6 ZI	grey	x	x		x	x	x																
	760663	Roto RFID 0,3 ZI	grey	x	x		x	x	x																
	760664	Roto RFID 2,5 Sintron	blue	x	x		x	x		x															
	760665	Roto RFID 1,0 Sintron	blue	x	x		x	x		x															
	760666	Roto RFID 0,6 Sintron	blue	x	x		x	x		x															
	760667	Roto RFID 0,3 Sintron	blue	x	x		x	x		x															
	760668	Roto RFID 2,5 PMMA	green	x	x		x	x			x	x		x											
	760669	Roto RFID 1,0 PMMA	green	x	x		x	x			x*	x	x		x										
	760670	Roto RFID 0,6 PMMA	green	x	x		x	x				x	x												
	760671	Roto RFID 0,3 PMMA	green	x	x		x	x			x*														
	760672	Roto RFID 2,0 Ti	black			x	x	x																	
	760673	Roto RFID 1,0 Ti	black			x	x	x																	
	760676	Diamond RFID 1,8	orange			x	x	x																	
	760677	Diamond RFID 1,4	orange			x	x	x																	
	760678	Diamond RFID 1,0	orange			x	x	x																	
	760679	Diamond RFID 0,4	orange			x	x	x																	
	760680	Roto RFID 2,5 SC	yellow		x		x	x						x		x									x*
	760681	Roto RFID 1,0 SC	yellow				x	x						x**											x*
	760682	Roto RFID 0,6 SC	yellow				x	x						x**											x*
	760683	Roto RFID 0,3 SC	yellow				x	x						x**											
	760684	Roto RFID 2,5 DMB DC	white	x	x	x	x	x																	
	760685	Roto RFID 1,0 DMB DC	white	x	x	x	x	x																	
	760686	Roto RFID 0,6 DMB DC	white	x	x	x	x	x																	
	760687	Roto RFID 0,3 DMB DC	white	x	x	x	x	x																	
	760688	Roto RFID 2,5 COCR	brown																						
	760689	Roto RFID 1,5 COCR	brown																						
	760690	Roto RFID 1,0 COCR	brown																						
	760692	Roto RFID 2,5 PMMA Denture	pink																						x
	760693	Roto RFID 1,5 PMMA Denture	pink																						x
	760694	Roto RFID 1,2 Drill	green																						
	760695	Roto RFID 1,0 Wax Denture	red																						x
	760696	Roto RFID 3,0 Wax Denture	red																						x
	181205	Roto RFID Calibration Pin	purple			x																			

* Processing only possible with Ceramill Matik and Ceramill Motion 2

** Processing only possible with Ceramill Matik

Fix Gel veneer holding gel

- Secure position.
- Easy removal.
- Burns clean and free.



Product Name	Quantity	Price
200-0000 FixGel	2 x 3.5ml	\$49

CeraFirm modeling fluid for ceramic

- Ceramic stays where it's put.
- Transformational layering accuracy.
- There are no limits.

Product Name	Quantity	Price
201-0000 CeraFirm Fluid Universal Ceramic Mixing Liquid	250ml	\$91
201-1000 CeraFirm Add Universal Ceramic Mixing Liquid	10ml	\$66
201-2000 CeraFirm Kit	1x10ml Add & 1x250ml Fluid	\$109



Super Peg ceramic firing support

- Precise firing fixation.
- Eliminates sandblasting entirely.
- Air blow or pluck out.



Product Name	Quantity	Price
8510 Super Peg II	4 x 12 cc syringe	\$47

Contour Check texture visual aid

- True texture visualization.
- Convenience in a liquid powder.
- Simply wipe or steam clean.



Product Name	Quantity	Price
7700 Contour Check, Silver	7ml	\$53

Clip simple abutment jig

- Precise pressure for bonding.
- Always accurate cementation.
- Universal for all implants.



Product Name	Quantity	Price
309-0000 Clip-Abutment Bonding Device	1	\$219

MS1 master die sealant

- Ultra-thin Sealant MS1 is absorbed right into the granular structure of the gypsum.
- Dries in seconds.

Product Name	Quantity	Price
7169 MS1 Master Die Sealer	1oz	\$39
7176 MS1 Master Die Sealer Refill w/empty bottle	8oz	\$69



MS3 master die separator

- Isolation of three materials: wax, ceramic and composite
- Ultra thin with no oily residue.
- Dries in seconds.



Product Name	Quantity	Price
7180 MS3 Master Die Separator	1oz	\$39
7183 MS3 Master Die Separator Refill w/ empty bottle	8oz	\$69



Model Sep gypsum isolation

- Easy gypsum separation.
- Clear liquid dries fast.
- Leaves no film residue.



Product Name	Quantity	Price
7065 Model-Sep Spray Bottle	8oz	\$58
7066 Model-Sep Refill w/ empty spray bottle	32oz	\$119

Surefit vibrant die spacer

- Accurate cement spacing.
- Excellent flow properties.
- Dries in seconds.



Product Name	Quantity	Price
7151 Surefit Die Spacer, Blue	0.5oz	\$33
7152 Surefit Die Spacer, Grey	0.5oz	\$33
7153 Surefit Die Spacer, Thinner	0.5oz	\$33
7156 Surefit Die Spacer, Dentine	0.5oz	\$33
7158 Surefit Die Spacer, Red	0.5oz	\$33
7159 Surefit Die Spacer, Gold	0.5oz	\$33
7160 Surefit Die Spacer, Silver	0.5oz	\$33
7162 Surefit Die Spacer Kit, Gold/ Silver/ Thinner	15ml	\$69
7163 Surefit Die Spacer KIT, Red/ Blue/ Thinner	15ml	\$69

Surfactant wetting agent/debubblizer

- Highly defined smooth casts.
- All gypsums and investments.
- Pressing wax or casting gypsums



Product Name	Quantity	Price
7060 Surfactant Trigger Spray	8oz	\$41
7061 Surfactant Refill w/ empty spray	32oz	\$84

True Prep prep shaded die spacer

- For All-Ceramic crowns and bridges.
- Shades master die.
- Available in 9 stump colors.



Product Name	Quantity	Price
7500 True Prep Kit, HD 1-9		\$401
7501 True Prep Spacer HD1	7ml	\$49
7502 True Prep Spacer HD2	7ml	\$49
7503 True Prep Spacer HD3	7ml	\$49
7504 True Prep Spacer HD4	7ml	\$49
7505 True Prep Spacer HD5	7ml	\$49
7506 True Prep Spacer HD6	7ml	\$49
7507 True Prep Spacer HD7	7ml	\$49
7508 True Prep Spacer HD8	7ml	\$49
7509 True Prep Spacer HD9	7ml	\$49
7510 True Prep Thinner	7ml	\$49

Dura Pattern light cure modeling composite

- Light cure efficiency.
- One-component.
- Flows like molten wax
- Cast patterns, attachments or implant jigs



Product Name	Quantity	Price
3000 DuraPattern LC Flow	3 x 3g Syringes	\$109

Versaglu dental cyanoacrylate

- Highly refined formula.
- Rapid structural bond.
- Unmatched tensile strength.
- Varying viscosities to suit the situation



Product Name	Quantity	Price
2025 Versaglu PX-4 Thin Adhesive	1oz	\$24
2026 Versaglu PX-40 Medium Adhesive	1oz	\$24
2027 Versaglu PX-400 Thick Adhesive	1oz	\$24

ZBase die blockout material

- Corrects Die anomalies.
- Full bodied and controlled.
- Flows like wax.
- Cure with ZBase accelerator



Product Name	Quantity	Price
2150 ZBase Blockout Material	10g	\$21
2151 ZBase Accelerator Spray	2oz	\$15
2152 ZBase Accelerator Refill w/ spray bottle	16oz	\$59

Speed Polish for zirconia and LiSi polish

- Polish ZrO₂ and LiSi₂ in seconds.
- Dust-free pickup.
- Made from scratch in small batches.

Recommended: #9 std stiff Abbott Robinson bristle brush



Product Name	Quantity	Price
8525 Speed Polish	27g	\$79

TriVest press plunger

- Disposable plungers ø13mm
- Rounded on both sides allows bi-directional pressing
- Enormous strength
- Ceramic universal



Product Name	Quantity	Price
222-0050 TriVest Plungers	13 mm PK/50	\$99

SureSep removable die spacer

- Spacer & Sep in one.
- Doubles as separator for acrylics, composites and metal abutments to stone
- Easily peels right off.
- Thixotropic coating – 20 microns



Product Name	Quantity	Price
7192 Sure Sep Composite Release	0.5oz	\$42

Instrum composite instruments

- Balanced handling
- Comfortable build-up
- Magnetic tip interchange
- Three tips:
 - > natural hair for build-up
 - > precision point to create structure and tattoo internal colour effects
 - > ball point to mimic gingiva and create an occlusal structure



Product Name	Quantity	Price
209-1000 Instrum Composite 1x Composite Build, 1x Composite Tatoo, 1x Composite Ball	Full Kit	\$199
209-1001 Instrum Composite Build	2/pk	\$91
209-1002 Instrum Composite Tatoo	2/pk	\$61
209-1003 Instrum Composite Ball	2/pk	\$73



Itsoclear clasps

Use Itsoclear clasps for: repairs, acrylic partials, cast partials with clear clasps in less than 3 minutes, replace metal clasps on existing partials, acrylic Nesbit with clear clasps.

- Crystal clear, will not lose clarity
- Non-abrasive
- Easily adjustable
- No injection machines or hand-held injectors required



Product Number	Description	Quantity	Price
1009151	Molar clasps 2.2mm thick	6 pack	\$74
1009150	Universal clasps 1.8mm thick	6 pack	\$74
1009153	Molar clasps 2.2 mm thick	36 pack	\$409
1009152	Universal clasps 1.8mm thick	36 pack	\$409

Tecnoflask microwave flask

- Acrylic material designed for use in microwave
- Cannot be damaged by exertion of pressure
- Polymerization time shortened to four minutes
- No isolating material needed between the muffle and the plaster
- Does not require a specific resin



Product Number	Description	Quantity	Price
1009199	Tecnoflask press board PL-100	each	\$139
1009200	Tecnoflask muffle	each	\$289
1009198	Tecnoflask kit	each	\$409

Clearmet clear thermoplastic polyamide

- Stain-resistant, odor-free material
- Monomer-free
- Clasps are adjustable with a warm three-prong plier, making it easy to adjust chairside
- Can be relined and repaired in saddle areas
- Clasps are tooth bearing, not tissue bearing
- 90% into undercut creates excellent retention



Product Number	Description	Quantity	Type*	Price
1009225	Small 25.5mm x 45mm	5 tubes/pack	S	\$82
1009227	Medium 25.5mm x 77mm	5 tubes/pack	S	\$99
1009229	Dental appliance	1kg	SO	\$399

*S = Stock
SO = Special Order

Versacryl



Versacryl Thermo-Soft Adjustable Reline

Product Number	Description	Quantity	*Type	Price
1014024	Sealer	0.5oz	S	\$30
1014025	Sealer	2oz	SO	\$109
1014021	Powder – pink	2oz	S	\$56
1014023	Powder – pink	1lb	SO	\$299
1014020	Powder – clear	2oz	S	\$56
1014022	Powder – clear	1lb	SO	\$299
1014010	Softening liquid	60ml	S	\$129
1014015	Softening liquid	1qt	SO	\$439

*S = Stock SO = Special Order

Versacryl Thermo-Soft Adjustable Reline

Self-cure

Product Number	Description	Quantity	*Type	Price
1014013	Self-cure pre-mixed	60ml	S	\$128
1014018	Self-cure pre-mixed	1qt	SO	\$439
1014011	Self-cure hardener	60ml	S	\$128
1014016	Self-cure hardener	1qt	SO	\$439
1014033	Self-cure reline kit	each	S	\$309
1014031	Self-cure lab kit	each	S	\$309

*S = Stock SO = Special Order

Versacryl Thermo-Soft Adjustable Reline

Heat-cure

Product Number	Description	Quantity	*Type	Price
1014014	Heat-cure pre-mixed	60ml	S	\$139
1014019	Heat-cure pre-mixed	1qt	SO	\$469
1014012	Heat-cure hardener	60ml	S	\$129
1014017	Heat-cure hardener	1qt	SO	\$439
1014034	Heat-cure reline kit	each	S	\$309
1014032	Heat-cure lab kit	each	S	\$309

*S = Stock SO = Special Order



DIAMOND D® High-Impact acrylic

Liquid

Product Number	Description	Quantity	*Type	Price
1013017	Self cure liquid	8oz (250 ml)	S	\$39
1013018	Self cure liquid	1qt (950 ml)	S	\$69
1013014	Heat cure liquid	8oz (250 ml)	S	\$31
1013015	Heat cure liquid	1qt (950 ml)	S	\$74
1013093	Diamond D 20-minute liquid	8oz (250 ml)	S	\$34
1013094	Diamond D 20-minute liquid	1qt (950 ml)	SO	\$84

*S = Stock SO = Special Order

Powder Heat-cure

Product Number	Description	Quantity	*Type	Price
1013020	Heat-cure powder Original	1lb	S	\$85
1013021	Heat-cure powder Original	5lb	S	\$339
1013022	Heat-cure powder Original	25lb	S	\$1179
1013024	Heat-cure powder Light	1lb	S	\$85
1013025	Heat-cure powder Light	5lb	SO	\$319
1013026	Heat-cure powder Light	25lb	SO	\$695
1013028	Heat-cure powder Light reddish pink	1lb	S	\$85
1013029	Heat-cure powder Light reddish pink	5lb	SO	\$319
1013030	Heat-cure powder Light reddish pink	25lb	SO	\$1095
1013032	Heat-cure powder Chroma Essence	1lb	SO	\$79
1013033	Heat-cure powder Chroma Essence	5lb	SO	\$319
1013034	Heat-cure powder Chroma Essence	25lb	SO	\$1095
1013105	Heat-cure powder Dark Veined	1lb	SO	\$79
1013106	Heat-cure powder Dark Veined	5lb	SO	\$319
1013107	Heat-cure powder Dark Veined	25lb	SO	\$1095

*S = Stock SO = Special Order

Powder and Liquid Self-cure

Product Number	Description	*Type	Price
1013073	1lb Powder & 8oz liquid - Original	S	\$112
1013074	1lb Powder & 8oz liquid - Light	SO	\$113
1013075	1lb Powder & 8oz liquid - Light reddish pink	SO	\$113
1013115	1lb Powder & 8oz liquid - Dark Veined	SO	\$113
1013077	5lbs Powder & 1qt liquid - Original	S	\$418
1013078	5lbs Powder & 1qt liquid - Light	S	\$418
1013079	5lbs Powder & 1qt liquid - Light reddish pink	S	\$418
1013116	5lbs Powder & 1qt liquid - Dark Veined	SO	\$399
1013081	25lbs Powder & 4qts Liquid - Original	SO	\$1279
1013082	25lbs Powder & 4qts Liquid - Light	SO	\$1279
1013083	25lbs Powder & 4qts Liquid - Light reddish pink	SO	\$1279

*S = Stock SO = Special Order

Powder Self-cure

Product Number	Description	Quantity	*Type	Price
1013050	Self-cure powder Original	1lb	S	\$99
1013051	Self-cure powder Original	5lb	SO	\$369
1013052	Self-cure powder Original	25lb	SO	\$1189
1013054	Self-cure powder Light	1lb	S	\$99
1013055	Self-cure powder Light	5lb	SO	\$369
1013056	Self-cure powder Light	25lb	SO	\$1189
1013058	Self-cure powder Light reddish pink	1lb	S	\$99
1013059	Self-cure powder Light reddish pink	5lb	SO	\$369
1013109	Self-cure powder Dark Veined	1lb	SO	\$89
1013110	Self-cure powder Dark Veined	5lb	SO	\$369
1013111	Self-cure powder Dark Veined	25lb	SO	\$1189

*S = Stock SO = Special Order

Powder and Liquid Heat-cure

Product Number	Description	*Type	Price
1013061	1lb Powder & 8oz liquid - Original	S	\$95
1013062	1lb Powder & 8oz liquid - Light	S	\$95
1013063	1lb Powder & 8oz liquid - Light reddish pink	S	\$95
1013064	1lb Powder & 8oz liquid - Chroma Essence	SO	\$89
1013112	1lb Powder & 8oz liquid - Dark Veined	SO	\$89
1013065	5lbs Powder & 1qt liquid - Original	S	\$369
1013066	5lbs Powder & 1qt liquid - Light	SO	\$329
1013067	5lbs Powder & 1qt liquid - Light reddish pink	SO	\$329
1013068	5lbs Powder & 1qt liquid - Chroma Essence	SO	\$329
1013113	5lbs Powder & 1qt liquid - Dark Veined	SO	\$329
1013069	25lbs Powder & 4qts liquid - Original	S	\$1199
1013070	25lbs Powder & 4qts liquid - Light	SO	\$1089
1013071	25lbs Powder & 4qts liquid - Light reddish pink	SO	\$1089
1013072	25lbs Powder & 4qts Liquid - Chroma Essence	SO	\$1089
1013114	25lbs Powder & 4qts Liquid - Dark Veined	SO	\$1089

*S = Stock SO = Special Order

Various Items

Product Number	Description	Quantity	*Type	Price
1013047	Mylar pouches	100/pack	SO	\$29
1013048	Mylar pouches	250/pack	S	\$56
1010170	Ultra Separator	1 gallon	S	\$69
1010175	Ultra Separator	1 quart	SO	\$28
1010180	Wax gloss	1 quart	SO	\$50
1013120	Ultrawash	1 gallon	SO	\$56

*S = Stock SO = Special Order



Enamelite Low-Fusing Fluorescent Ceramic Spray Glaze

- Guarantees shade-matching under natural and artificial light for zirconia, lithium-disilicate, and lithium-silicate restorations
- Stain and glaze in one firing
- Won't puddle at margins
- Glaze multiple thin, uniform, & accurate coats IN SECONDS
- Faster drying than brush-on applications
- Compatible with ALL dental ceramics



Product Name	Quantity	Price
5100181	46g/ 1.6oz	\$109

Enamelite CAD/CAM Scanning Spray



Product Name	Quantity	Price
5100176	88g/ 3.1oz	\$109

Mizzy Heatless Wheels - RED

- Aluminum oxide with ceramic coating
- Use on non and semi precious alloys, porcelain and zirconia
- No need to switch back to silicone carbide wheel (Mizzy grey) for final finishing
- 50/box \$45



Product Number	Size	Dimension	*Type
6131608	#8	3.1 x 19mm	S
6131605	#5	3.1 x 22.2mm	S
6131603	#3	2.4 x 25.4mm	S
6131602	#2	3.1 x 25.4mm	S
6131600	#1	4.8 x 25.4mm	S

*S = Stock SO = Special Order

Mizzy Heatless Wheels - Grey

- Silicone carbide
- Rapid reduction through final shaping
- Use on non-precious, semi-precious and precious alloys, porcelain, acrylic, and denture teeth
- Use up to 20,000 RPM
- 50/box \$42



Product Number	Size	Dimension	*Type
6131500	#15	2.4 x 12.7mm	SO
6131400	#14	3.1 x 12.7mm	SO
6131300	#13	4.8 x 12.7mm	SO
6131200	#12	2.4 x 15.9mm	SO
6131100	#11	3.1 x 15.9mm	SO
6131000	#10	4.8 x 15.9mm	SO
6130900	#9	2.4 x 19mm	SO
6130800	#8	3.1 x 19mm	SO
6130700	#7	4.8 x 19mm	SO
6130600	#6	2.4 x 22.2mm	SO
6130500	#5	3.1 x 22.2mm	SO
6130400	#4	4.8 x 22.2mm	SO
6130300	#3	2.4 x 25.4mm	SO
6130200	#2	3.1 x 25.4mm	SO
6130100	#1	4.8 x 25.4mm	SO

*S = Stock
SO = Special Order

Mizzy Heatless Wheels - White

- Aluminum oxide
- Use on porcelain or semi-precious and precious alloys
- 50/box \$45



Product Number	Size	Dimension	*Type
6151500	#15	2.4 x 12.7mm	SO
6151400	#14	3.1 x 12.7mm	SO
6151300	#13	4.8 x 12.7mm	SO
6151200	#12	2.4 x 15.9mm	SO
6151100	#11	3.1 x 15.9mm	SO
6151000	#10	4.8 x 15.9mm	SO
6150900	#9	2.4 x 19mm	SO
6150800	#8	3.1 x 19mm	SO
6150700	#7	4.8 x 19mm	S
6150600	#6	2.4 x 22.2mm	S
6150500	#5	3.1 x 22.2mm	S
6150400	#4	4.8 x 22.2mm	S
6150300	#3	2.4 x 25.4mm	S
6150200	#2	3.1 x 25.4mm	S
6150100	#1	4.8 x 25.4mm	S

*S = Stock SO = Special Order

PRO-FORM

Proform sheets are available in a large variety of sizes and types. If an item you are looking for is not listed in this section please let us know and we can place a special order for you.

Mouthguard Laminates

(5" x 5" square - 4mm thick)
12/box **\$114**





*S = Stock
SO = Special Order

Product Number	Description	*Type
9596920	Clear	S
9598060	Blue	S
9598710	Two-color Assorted	S
9598730	Black	S
9600170	Three-color assorted	S
9600180	One color- assorted	S
9600181	Fun Guards – Assorted	S
9600800	Green	S
9600820	Maroon	SO
9600850	Orange	S
9600870	Purple	SO
9600886	Pink	SO
9600889	Breast Cancer pink	SO
9600910	Red	SO
9601020	White	SO
9601040	Yellow	S



Proform Splint




Product Number	Description	Quantity	*Type	Price
9614780	.060 (1.5mm)	50/pack	SO	\$50
9614790	.080 (2.0mm)	50/pack	SO	\$64
9614800	.020 (0.5mm)	50/pack	S	\$35
9614810	.030 (0.8mm) temp	50/pack	SO	\$38
9614816	.030 (0.8mm) clear retainer	50/pack	SO	\$145
9614820	.100 (2.5mm)	50/pack	SO	\$73
9614830	.040 (1.0mm)	50/pack	SO	\$44
9614960	.060 (1.5mm)	25/pack	S	\$29
9614980	.080 (2.0mm)	25/pack	S	\$35
9615000	.020 (0.5mm)	25/pack	SO	\$17
9615020	.030 (0.8mm) temp	25/pack	SO	\$35
9615025	.030 (0.8mm) clear retainer	25/pack	SO	\$85
9615040	.100 (2.5mm)	25/pack	SO	\$45
9615060	.040 splint (1.0 mm)	25/pack	SO	\$24

*S = Stock SO = Special Order

Proform	Product Number	Description	Quantity	*Type	Price
	9616210	Base Material .060 (1.5mm)	25/pack	SO	\$23
	9616230	.080 (2.0mm)	25/pack	SO	\$28
	9617390	Tray .125 (3.2 mm)	25/pack	SO	\$38
	9617410	.150 (3.8mm)	25/pack	SO	\$42

*S = Stock
SO = Special Order

Proform	Product Number	Description	Quantity	*Type	Price
	9602550	Niteguard	12/box	S	\$107
	9603580	Dual Laminates – E-Gasket Clear 3mm	12/box	S	\$109
	9604750	Pink	12/box	SO	\$109

Proform	Product Number	Description	Quantity	*Type	Price
	9596980	Soft EVA .040 (1.0mm)	25/box	SO	\$32
	9597090	.060 (1.5mm)	25/box	SO	\$41
	9597150	.080 (2.0mm)	25/box	SO	\$44
	9611240	.120 (3mm)	25/box	S	\$43
	9613590	Coping .020 (0.5mm)	25/box	SO	\$16
	9613610	.030 (0.8mm)	25/box	SO	\$17
	9613630	.040 (1.0mm)	25/box	SO	\$24
	9613720	Crown & Bridge .020 (0.5mm)	25/box	SO	\$25

Tray Always Disposable Impression Trays – Edentulous



*S = Stock
S0 = Special Order

Product Number	Size	Description	Quantity	*Type	Price
0921599	#2	lower small	12/bag	S0	\$41
0921600	#2	upper small	12/bag	S0	\$41
0921700	#3	upper medium	12/bag	S	\$41
0921701	#3	lower medium	12/bag	S	\$41
0921702	#5	upper large	12/bag	S	\$41
0921703	#5	lower large	12/bag	S	\$41
0921704		variety bag – 2 each size	16/bag	S	\$41
0921706	#7	upper extra large	12/bag	S	\$41
0921707	#7	lower extra large	12/bag	S	\$41

Lathe Burs

Keystone Spiral Cut Carbide Lathe Burs

- Require less pressure and cut smoothly
- Minimum heat generation
- Hand ground for optimal quality



Product Name	Description	Price
1201680	A 1/2	\$39
1201690	B 1/2	\$39
1201700	C 1/2	\$39
1201710	D 1/2	\$39
1201720	E 1/2	\$39

Keystone Diamond Cut Carbide Lathe Burs - Coarse Maxi Cut

- Coarse maxi cut for maximum reduction
- Use on acrylic and metals



Product Name	Description	Price
1202181	A 1/2	\$48
1202121	B 1/2	\$48
1202123	D 1/2	\$42
1202125	A 3/8	\$39
1202128	D 3/8	\$39

Keystone Diamond Cut Carbide Lathe Burs - Fine Mini Cut


- Fine mini cut for minimum reduction
- Smooth finish on acrylic, composite, porcelain and metal



Product Name	Description	Price
1202154	D 1/2	\$59
1202159	D 3/8	\$39

Keygold Carbide Lathe Burs

- Titanium nitride coating
- Require less pressure and cut smoothly
- Minimum heat generation
- Hand ground for optimal quality



Product Name	Description	Price
1202230	A 1/2	\$48
1202235	B 1/2	\$48
1202245	D 1/2	\$48
1202255	A 3/8	\$39
1202270	D 3/8	\$39

Hatho Miniature Handpiece Scotch Brite Brushes

5,000 RPM
12/pack \$34

Product Name	Description
1670016	Coarse
1670092	Medium
1670091	Fine



Hatho Scotch Brite Polishing Lathe Discs

- Ideal for softer materials like denture acrylic and denture reline material
- Wood core encased in plastic allows for easy mounting/dismounting



Product Name	Description	Price
FINE		
1670012	4" dia., 1" thick, 2,800 rpm	\$27
1670013	3" dia., 1" thick, 1,400 rpm	\$22
1670014	4" dia., 1/2" thick, 2,800 rpm	\$20
1670015	3" dia., 1/2" thick, 1,400 rpm	\$16

Product Name	Description	Price
MEDIUM		
1670008	4" dia., 1" thick, 2,800 rpm	\$29
1670009	3" dia., 1" thick, 1,400 rpm	\$23
1670010	4" dia., 1/2" thick, 2,800 rpm	\$19
1670011	3" dia., 1/2" thick, 1,400 rpm	\$17

Product Name	Description	Price
COARSE		
1670004	4" dia., 1/2" thick, 2,800 rpm	\$27
1670005	3" dia., 1/2" thick, 1,400 rpm	\$18
1670006	4" dia., 1/2" thick, 2,800 rpm	\$27
1670007	3" dia., 1/2" thick, 1,400 rpm	\$18

DETAX

FREEFORM[®] Light-cure Modeling Resin



- For the fabrication of surgical and x-ray templates, splints, base plates and transfer jigs
- Radio transparent
- Heat sterilizable
- Clear-transparent

Product Number and Name	Unit	Price
02204 Starter set plast fixgel bond coat	60g jar 3g syringe 5ml bottle 10ml bottle	\$149
02542 plast	60g jar	\$84
02829 fixgel	3g syringe	\$26
02702 bond	5ml bottle	\$16
02549 coat	10ml bottle	\$26

easyform light-cure modeling material



- Methacrylate free
- Available in gel or paste
- Burns out without residue
- Can be combined with regular or light-cured wax
- Translucent red

Product Number and Name	Unit	Price
03376 Gel application needles	3 x 3g syringes 10pcs	\$75
03375 Paste	30g jar	\$139

FLEXSTONE[®] Plus insulating and modeling silicone



- A-silicone
- 85 Shore A final hardness
- Recommended by Ivoclar-Vivadent for SR-Ivocap system
- Base – blue, Catalyst – white (1:1)

Product Number and Name	Unit	Price
02383 base + catalyst with stand mixing bowl	2 x 160ml bottles 1pc.	\$139
02329 base + catalyst mixing bowls	2x 1000ml 4pcs.	\$649

blue eco lab putty



- Heat-resistant up to 200C
- Final hardness 86 Shore A.
- Colour: Base blue, Catalyst beige

Product Number and Name	Unit	Price
02038 Standard pack base catalyst measuring scoops	800g jar 800g jar 2 pcs.	\$62
02467 Mega pack base catalyst measuring scoops	5kg tub 5kg tub 2 pcs.	\$379

blue eco scan® scannable lab putty



- Final hardness 85 Shore A.
- Colour: Base light grey, Catalyst white
- Eliminates need for CAD/CAM spray

Product Number and Name	Unit	Price
03118 Standard pack base catalyst measuring scoops	800g jar 800g jar 2 pcs.	\$61

blue eco stone lab putty



- Heat-resistant up to 200C
- Final hardness to 96 Shore A
- Colour: Base dark blue, Catalyst white

Product Number and Name	Unit	Price
03117 Standard pack base catalyst measuring scoops	800g jar 800g jar 2 pcs.	\$61
02517 ECO pack base catalyst measuring scoops	1400g jar 1400g jar 2 pcs.	\$94
02985 Mega pack base catalyst measuring scoops	5kg tub 5kg tub 2 pcs.	\$359

Implant Mask rigid gingival mask for implants



- Available in scan version
- Final hardness 70 Shore A
- Colour: gingiva

Product Number and Name	Unit	Price
02529 Standard pack mini-mix 1:1 mixing cannulas, light green separating liquid	2 x 10ml cartridges 6 pcs. 15ml bottle	\$42
03374 Standard pack scan mini-mix 1:1 mixing cannulas, light green separating liquid	2 x 10ml cartridges 6 pcs. 15ml bottle	\$46
03552 Cartridge pack scan cartridges mixing cannulas, blue separating liquid	2 x 50ml 10 pcs. 15ml bottle	\$72



Esthetic Mask flexible gingival mask



- Final hardness 45 Shore A
- Colour: gingiva

Product Number and Name	Unit	Price
03113 Cartridge packing cartridges mixing cannulas, yellow separating liquid	2 x 50ml 10 pcs. 15ml bottle	\$79

PERMALINER silicone soft denture liner



- Permanent soft reline material
- Heat cure
- Bonds to acrylics without bonding material
- Methacrylate free

Product Number and Name	Unit	Price
02092 Standard pack	30g jar	\$79
02146 Large pack	70g jar.	\$169

Molloplast® B silicone reline



- Lab processed permanent soft reline material
- Heat cure

Product Number and Name	Unit	Price
03000 Standard pack	45g jar	\$139
03001 Large pack	170g jar	\$299
03002 Combi packing jar	45g	\$159
Primo Adhesive		
03004 Primo adhesive	15ml bottle	\$36

mollosil® plus silicone soft reline cold cure



- Long-term soft relining material
- Chairside or indirect use
- Primer with colour indicator for ease of use

Product Number and Name	Unit	Price
02353 Standard pack Automix2 mixing cannulas, pink	50ml cartridge 7 pcs.	\$346
mollosil plus Primer	5ml bottle	
mollosil plus Polish	2 x 7ml bottles	
grinding bur, pointed	1 pc.	
grinding bur holder, pointed	1 pc.	
Pipettes	2 pc.	
02354 Refill pack Automix2	50ml cartridge	\$169
02440 mollosil plus Primer	5ml bottle	\$39
02441 mollosil plus Polish	2 x 7ml bottles	\$49



Snow Scan Brush-On CAD/CAM Powder

- Does not clump like spray powders
 - Thin layers = 10x more accuracy
 - No overspray = 30x more scan sessions
- 10g - \$199**



	Snow Scan Powder 10g	Scan Spray Type 72ml
Accuracy	1µ	12µ
Usage	3000 abutments	50-60 abutments

NEW! Detaseal Hydroflow VPS - Precision Impression and Bite Registration materials



Light and extra light Snap set

2 x 50ml with six yellow tips


Regular Set:

- Working time: 2:00
- Setting time: 2:00

Fast Set:

- Working time: 1:00
- Setting time: 1:30

\$69



Mono #02754

2 x 50ml with six green tips

Set:

- Working time: 1:30
- Setting time: 2:00

\$69



Mono - Jumbo Pack #03373

1 x 380ml with ten mixers

Set:

- Working time: 1:30
- Setting time: 2:00

\$119



Putty #02727

2 x 250ml

Set:

- Working time: 1:30
- Setting time: 2:00

\$72



Heavy - Jumbo Pack #02729

1 x 380ml with ten mixers

Set:

- Working time: 1:30
- Setting time: 2:00

\$119



Soft2 Putty with mango aroma #02880

2 x 500ml

Set:

- Working time: 1:15
- Setting time: 1:30

\$119



Detaseal Bite 90 Shore A #02139

2 x 50ml with six rose tips

Set:

- Working time: 0:45
- Setting time: 1:30

\$72



NEW! Detaseal Hydroflow VPS - Precision Impression and Bite Registration materials



Greenbite Colour with colour shift and green apple flavour
45 Shore A #02207

2 x 50ml with six green tips

Set:

- Working time: 0:30
- Setting time: 0:50

\$72



Milkbite scannable with ice cream aroma
90 Shore A #03380

2 x 50ml with six rose tips

Set:

- Working time: 0:30
- Setting time: 0:50

\$69



Greenbite Apple with green apple flavour
45 Shore A #02641

2 x 50ml with six green tips

Set:

- Working time: 0:30
- Setting time: 1:00

\$72



Glassbite clear
80 Shore A #02571

1 x 50ml with four green tips

Set:

- Working time: 0:30
- Setting time: 1:00

\$39





PRODUCER OF SPECIAL ALLOYS SINCE 1975

Magnum HBA Chrome-Cobalt for partial frameworks

- Composition: 64% chrome, 29% cobalt, 6% Molybdenum, others (C, Si, Fe, Mn)
- Density: 8.4g/cm³
- Melting point: 1460C
- Yield load strength: 580 MPa
- CTE @500C: 14.1
- Part # 40160030A

1kg - \$289



Tilite-P Non Precious alloy for C&B frameworks

- CTE @500C: 14.0
- CTE @600C: 14.1
- Melting range: 1204-1302C
- Casting temp: 1329C



5oz - \$199



Abbott- Robinson brushes

Straight HP mounted bristle brushes



Size	Diameter
No.9	9/16"
No.11	11/16"
No. 12	3/4"

Brush Style	Bristle Colour
Soft	Brown
Standard-Stiff	Sliver-Grey
Extra-Stiff	Black

Product Number	Description	Quantity
06044	#9 X-Stiff	144
06014	#9 Std-Stiff	144
05994	#9 Soft	144
06104	#11 X-Stiff	144
06048	#11 Std-Stiff	144
06064	#11 Soft	144
06194	#12 X-Stiff	144
06164	#12 Std-Stiff	144
06144	#12 Soft	144
		Packs of 144=\$158

Product Number	Description	Quantity
06041	#9 X-Stiff	12
06011	#9 Std-Stiff	12
05991	#9 Soft	12
06101	#11 X-Stiff	12
06081	#11 Std-Stiff	12
06061	#11 Soft	12
06191	#12 X-Stiff	12
06161	#12 Std-Stiff	12
06141	#12 Soft	12
		Packs of 12=\$22

ALL ARTEX ARTICULATORS AT ONE SIGHT

Technical Data Chart	Non Arcon-Design			Arcon-Design	
	Artex BN	Artex CN	Artex CT	Artex CPR	Artex CR
SCI (Sagittal condylar inclination)	35°	35°	-15° to +60°	-20° to +60°	-20° to +60°
Bennett-angle (HCI)	15°	0° to +20°	0° to +20°	-5° to +30°	-5° to +30°
Protrusion	-	-	-	-	0 to 6mm
Retrusion	-	-	-	0 to 2mm	0 to 2mm
ISS (Immediate Sideshift)	-	-	-	-	0 to 1.5mm (per side)
Distraction	-	-	-	0 to 3mm	0 to 3mm
Centric design	Centric click	Centric click	Centric click	Centric lever for semi-axes	Centric lever for semi-axes
Upper and lower articulator arms prevented from coming apart inadvertently in open centric by:	Centric plate	Centric plate	Centric plate	Arcon-Clip	Arcon-Clip
Arbitrary pins for direct transfer with Artex facebow	no	no	yes	yes	yes

Artex Facebow

- Cranium/axis relation determined in corresponding to the patient's anatomy in just two minutes
- Quick and secure the 3-D universal joint
- Cranium-axis oriented transfer of maxillary model to the articulator transfers the real jaw position precisely to the articulator and consequently reduces grinding time at the chairside
- Reliable reproducibility of arbitrary axis position using cushioned Leipzig Nasion
- Product #218600 **\$1235**



Giroform Pin Drill

- Laser beam for easy drill positioning
- Plate holder smoothly adjustable thus ensuring safe and fast operation
- Plate holder is secured magnetically and automatically when drilling starts
- Automatic drill advance guarantees identical boreholes in the plates (0.5 seconds per drill cycle)
- Ergonomically positioned tip-on buttons for a high degree of operating convenience and easy activation of drilling
- Guide grooves ensure easy drill change
- Universal plate holder for all sizes and shapes of Giroform baseplates with anti-rotation protection
- Product #176701 **\$3450**



Smartbox X2 Plaster Dispenser

- Ensures high degree of dosage precision resulting in more accurate models
- Improves the consistency of the water/stone mixture

Technical data:

- Dimensions: 335 x 240 x 590mm
- Weight: 16.5kg
- Electrical connected value: 100/115/230 V, 50/60Hz
- Power Output: 95W
- Capacity Powder: appr. 8kg
- Time to dispense 20ml: 100g appr. 30s
- Product #116170 **\$4450**



Smartmix X2 Vacuum Mixer

- Multi-functional dial / button for quick and easy program selection and programming
- Homogeneous and consistent, reproducible mixtures with optimum parameters
- Quick access to all mixing parameter: Mixing duration, speed, direction of rotation, interval times, pre-mixing, pre-vacuum, post-vacuum and program name
- Space saving, can be used as a free-standing or wall unit

Technical data:

- Dimensions: 250 x 160 x 350mm (without stand)
- Weight: 8.5kg
- Electrical connected value: 100/115/230V, 50/60Hz
- Power Output: 210W
- Delivery rate: 15.8l/min
- Vacuum power: -800mbar
- Product #1157000 **\$2795**



Noflame Plus

The "electric Bunsen burner"

- Mobile and immediately operational, no gas required
- Handling as usual, no adjustment required
- Eliminates the risks of getting burned or causing fire
- No heat emission into the environment; does not consume oxygen
- Clean, i.e. no soot particles on instrument and/or material
- Saves 75% of energy expenses; saves the environment
- Ideal conditions, to be used even at the dental office
- Product #116250V115 **\$699**



Smartwax Duo

- Duo-option, i.e. two different handpieces are ready for use at the same time
- 3 individually programmable temperatures for each handpiece
- Simple and safe exchange mechanism of the modeling tips
- Fast operational and heated instruments
- Anti-twist cables that are insensitive to heat
- Gold plated copper tips for better heat conductivity

Technical data:

- Dimensions: 130 x 150 x 50mm
- Electrical connected value: 110-230V 50/60Hz
- Power Output: 12W
- Length of the handpiece cord: 1.8m
- Output voltage: 6V
- Temp.-range: 50-220°C/122-428°F - switchable
- Product #116270V115 **\$635**



Procure 300 Light Oven

The ProCure 300 Light Oven is an advanced third generation light curing unit that provides superior lighting intensity and excellent depth of cure.

Product #
9634590

Price
\$918



Pro-form® Machine III vacuum former

The Machine III has been a staple in Keystone Industries' long-lasting line of highly-recommended vacuum former machines. The price, features and reliability make it the first choice in dental vacuum formers. The machine has more power than any other dental vacuum former, which results in a more accurate and definite dental thermoplastic appliance.

Product #
7000330

Machine III – 110v

Price
\$699



Pro-form® Single Chambered Vacuum Former

The Single Chambered Vacuum Former is an economical answer that can still excellently fabricate most thermo-plastic appliances.

Product #
9594470

110v

Price
\$599 (special order)



recommended applications

Product Line	Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium	Peek	Page
PANTHER	■		■							40
WILDCAT		■								42
MARGAY		■								43
SERVAL	■		○					○		44
JUNGLE CAT	■		■			○		○		45
BOBCAT	○		■			○		○		46
JAGUAR			■							48
LEOPARD									■	49
DESERT CAT				■						50
LYNX			○	○	○	○	■	○		51
COUGAR			○	○	○		■			52
OCELOT					■					53
PUMA					■	○	○	○		53
CHEETAH					■					54
TIGER						■		○		55
LION C+M						○	■	○		56
LION F+XF						■	■			57

■ recommended ○ suitable

LAB SUPPLIES (ROTARY)

PANTHER

For Zirconium Oxide / Ceramic












Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
■		■					

- ☑ **Ultra-Premium Line for all-ceramic restorations**
 - ☑ High-gloss polish in just two steps
 - ☑ Perfect surface for antagonist protection
 - ☑ Long service life
- ☑ Roughened surfaces for immediate use
 - ☑ Various shapes for all situations
 - ☑ Soft bonding for smooth grinding



PANTHER M	Part # PANTHER100	PANTHER102	PANTHER104	PANTHER108
Grit 	medium			
	3 pack \$42.95	Single \$70.95	3 pack \$104.95	3 pack \$104.95
				
Sizes [mm]	Ø 5.5 x 18	Ø 26 x 2	Ø 17 x 2.4	Ø 15 x 1
Shank	104-HP	104-HP	104-HP	104-HP
○ optimal / maximum rpm	12,000/20,000	8,000/20,000	8,000/20,000	8,000/20,000


PANTHER XF	Part # PANTHER101	PANTHER103	PANTHER105	PANTHER109
Grit 	extra fine			
	3 pack \$42.95	Single \$70.95	3 pack \$104.95	3 pack \$104.95
				
Sizes [mm]	Ø 5.5 x 18	Ø 26 x 2	Ø 17 x 2.4	Ø 15 x 1
Shank	104-HP	104-HP	104-HP	104-HP
○ optimal / maximum rpm	12,000/20,000	8,000/20,000	8,000/20,000	8,000/20,000

Equipment	Part # PANTHER110	PANTHER106	PANTHER107	PANTHER111
	Single 139.95	10 pack \$26.95	10 pack \$26.95	
				
Sizes [mm]		Ø 17	Ø 19	
Shank	Diamond Polishing Paste (20 grams)	104-HP	104-HP	Panther edition kit
○ optimal / maximum rpm		12,000/15,000	12,000/15,000	Kit 475.95

Panther Abutment Polishing System

The surface of the implant abutment in the transmucosal area plays an important role. If the surface is too rough, the risk of increased plaque accumulation occurs. If the surface is too smooth, the fibroplasts of the peri-implant mucosa cannot "fuse/heal" optimally. Various studies have determined a threshold exists when plaque adhesion is low and fibroplasts can fuse/heal at the same time.

The ideal surface roughness is between (μm : Ra = 0.21 - 0.40).



PANTHER

For Zirconium Oxide / Ceramic



The PANTHER ABUTMENT SURFACE protocol can be used for titanium, ceramic or hybrid abutments by following the two steps below.

1. In the first step, excess glue and other materials should be removed with the PANTHER lense rough. This instrument has a coarse diamond grain.



2. In the second step, the desired residual roughness will be applied. The PANTHER lense smooth achieves an average roughness of approx. 0.34 μm .



PANTHER R		Part # PANTHER113	PANTHER115
Grit	Rough	Single \$79.95	3 pack \$104.95
Sizes [mm]		Ø 26 x 2	Ø 15 x 1
Shank		104-HP	104-HP
○ optimal / maximum rpm		8,000/20,000	8,000/20,000

PANTHER S		Part # PANTHER102	PANTHER108
Grit	Smooth	Single \$70.95	3 pack \$104.95
Sizes [mm]		Ø 26 x 2	Ø 15 x 1
Shank		104-HP	104-HP
○ optimal / maximum rpm		8,000/20,000	8,000/20,000

NEW			
PANTHER Spirals		Part # Panther 120	Panther 121
		Panther 120 - \$39.95	
		Panther 121 - \$95.95	
Grit	Medium	Grit	Fine
Sizes [mm]		Ø 20,3	Ø 20,3
○ optimal / maximum rpm		8,000/20,000	8,000/20,000



Panther Abutment Polishing Kit

Part # PANTHER112

Kit 265.95

LAB SUPPLIES (ROTARY)

WILDCAT

For Zirconium Oxide (unsintered)



Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
	■						

- Specially designed for green state zirconia applications
 - High durability
 - Smooth handling

WILDCAT | M

Grit medium

| Sizes [mm]

| Shank

○ optimal / maximum rpm

Part # **WILDCAT100** **WILDCAT102** **WILDCAT104**

10 pack \$37.95 **10 pack \$37.95** **10 pack \$42.95**



Ø 5.6 x 13

Ø 7 x 13

Ø 15 x 0.3

104-HP

104-HP

104-HP

6,000 / 20,000

6,000 / 20,000

6,000 / 20,000

WILDCAT | M/F

Grit medium/fine

| Sizes [mm]

| Shank

○ optimal / maximum rpm

Part # **WILDCAT106**

Single \$14.95

10 pack \$103.95



Ø 7 x 19

104-HP

8,000 / 20,000



Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
	■						

MARGAY

For Zirconium Oxide (unsintered)



- Specially designed for removal of milling sprues**
- Developed with dental technicians



MARGAY | M/F

Grit medium/fine

| Sizes [mm]

| Shank

⊖ optimal / maximum rpm

Part # **MARGAY100**

Single **\$15.95**

10 pack **\$139.95**

Ø 10 x 22.7

104-HP

8,000 / 20,000

The unique Weissling polisher was developed in close collaboration with experts and laboratories working directly with such practical applications.

The special features include the precise form of the polisher, which is perfectly tailored to the task and its two assembled abrasive material zones. The grey zone reduces zirconia quickly, and the pigment-free white zone is used for precise fine-polishing, with no risk of surface colouring or damage to the precious working part.

The result: precise, efficient and reduced risk of damage.



Grey zone: quick reduction of the milling sprues.



White Zone: Step two – same polisher. Risk-free contouring through reduced removal.



Two in One: fast, accurate and reliable processing.

LAB SUPPLIES (ROTARY)

SERVAL

For Zirconium Oxide (sintered)



Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
■		○					○

- ✓ Premium Line For Sintered Zirconia restorations
- ✓ Excellent wear and tear characteristics
- ✓ High concentration of diamond grit



SERVAL XC	Part #	SERVAL116	SERVAL146	SERVAL108	SERVAL126	SERVAL138
Grit		extra coarse				
		Single \$59.95	Single \$38.95	Single \$20.95	Single \$20.95	Single \$17.95
		10 pack \$149.95	10 pack \$149.95	10 pack \$149.95	10 pack \$149.95	10 pack \$139.95
						
Sizes [mm]		Ø 26 x 2	Ø 22 x 3.2	Ø 5.5 x 16	Ø 11.5 x 2.4	Ø 6 x 15
Shank		104-HP	104-HP	104-HP	104-HP	104-HP
○ optimal / maximum rpm		6,000 / 20,000	6,000 / 20,000	15,000 / 20,000	6,000 / 20,000	15,000 / 20,000

SERVAL M	Part #	SERVAL114	SERVAL144	SERVAL104	SERVAL122	SERVAL134
Grit		medium				
		Single \$59.95	Single \$38.95	Single \$20.95	Single \$20.95	Single \$17.95
		10 pack \$149.95	10 pack \$149.95	10 pack \$149.95	10 pack \$149.95	10 pack \$139.95
						
Sizes [mm]		Ø 26 x 2	Ø 22 x 3.2	Ø 5.5 x 16	Ø 11.5 x 2.4	Ø 6 x 15
Shank		104-HP	104-HP	104-HP	104-HP	104-HP
○ optimal / maximum rpm		6,000 / 20,000	6,000 / 20,000	15,000 / 20,000	6,000 / 20,000	15,000 / 20,000

SERVAL F	Part #	SERVAL112	SERVAL142	SERVAL100	SERVAL118	SERVAL130
Grit		fine				
		Single \$59.95	Single \$38.95	Single \$20.95	Single \$20.95	Single \$17.95
		10 pack \$149.95	10 pack \$149.95	10 pack \$149.95	10 pack \$149.95	10 pack \$139.95
						
Sizes [mm]		Ø 26 x 2	Ø 22 x 3.2	Ø 5.5 x 16	Ø 11.5 x 2.4	Ø 6 x 15
Shank		104-HP	104-HP	104-HP	104-HP	104-HP
○ optimal / maximum rpm		6,000 / 20,000	6,000 / 20,000	15,000 / 20,000	6,000 / 20,000	15,000 / 20,000

Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
■		■			○		○

JUNGLE CAT

For Ceramic



- ☑ **Our Premium Line**
- ☑ Ceramic bonded grinding stone for ceramic and Zirconia
- ☑ Special formulation for an extremely long service life
- ☑ Minimal residual roughness
- ☑ Steady grinding with low heat



JUNGLE CAT | C

Grit Coarse



Sizes [mm]
Shank
○ optimal / maximum rpm

Part #	JUNGLECAT100	JUNGLECAT106	JUNGLECAT114	JUNGLECAT132
	Single \$37.95	Single \$37.95	Single \$46.95	Single \$57.95
	Ø 5 x 13	Ø 4 x 10.5	Ø 15 x 3.5	Ø 22 x 2
	104-HP	104-HP	104-HP	104-HP
	15,000 / 20,000	15,000 / 20,000	10,000/20,000	10,000/20,000

JUNGLE CAT | M

Grit medium



Sizes [mm]
Shank
○ optimal / maximum rpm

Part #	JUNGLECAT118	JUNGLECAT102	JUNGLECAT104	JUNGLECAT108	JUNGLECAT110
	Single \$46.95	Single \$37.95	Single \$37.95	Single \$37.95	Single \$41.95
	Ø 12 x 6.5	Ø 5 x 13	Ø 6.5 x 13	Ø 4 x 10.5	Ø 12 x 2.4
	104-HP	104-HP	104-HP	104-HP	104-HP
	10,000 / 20,000	15,000 / 20,000	15,000 / 20,000	15,000 / 20,000	10,000 / 20,000

Part #	JUNGLECAT112	JUNGLECAT128	JUNGLECAT134
	Single \$46.95	Single \$59.95	Single \$57.95
	Ø 15 x 3.5	Ø 22 x 2.1	Ø 22 x 2
	900 / 104-HP	900 / 104-HP	900 / 104-HP
	10,000 / 20,000	10,000 / 20,000	10,000 / 20,000
	A / B	A / B	A / B

| matching mandrels please see page 54

Sizes [mm]
Shank
○ optimal / maximum rpm
Mandrels

LAB SUPPLIES (ROTARY)

BOBCAT





For Ceramic




Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
○		■			○		○

- ✓ **Our Value Line**
- ✓ Suitable for all kind of ceramics
- ✓ High durability
- ✓ Effective and quick polishing
- ✓ Proven diamond technology



BOBCAT C	BOBCAT124	BOBCAT160	BOBCAT152
Grit	coarse		
	Single \$53.95	Single \$28.95	Single \$34.95
			
Sizes [mm]	Ø 26 x 2	Ø 17 x 2.6	Ø 22 x 3.2
Shank	104-HP		
○ optimal / maximum rpm	6,000 / 20,000		

BOBCAT C	10 pack \$149.95	10 pack \$99.95	BOBCAT176
Grit	coarse		
	Single \$16.95	Single \$13.95	Single \$36.95
			
Sizes [mm]	Ø 5.5 x 18	Ø 11.5 x 2.4	Ø 22 x 3.5
Shank	104-HP		
○ optimal / maximum rpm	6,000 / 20,000		



Polisher surface grit C (coarse)



Polisher surface grit M (medium)



Polisher surface grit F (fine)



Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
○		■			○		○

BOBCAT


For Ceramic



BOBCAT M	BOBCAT128	BOBCAT168	BOBCAT156
Grit medium	Single \$53.95	Single \$28.95	Single \$34.95
			
Sizes [mm]	Ø 26 x 2	Ø 17 x 2.6	Ø 22 x 3.2
Shank	104-HP	104-HP	104-HP
○ optimal / maximum rpm	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

BOBCAT M	BOBCAT108	BOBCAT136	BOBCAT180
Grit medium	10 pack \$149.95 Single \$16.95	10 pack \$99.95 Single \$13.95	Single \$36.95
			
Sizes [mm]	Ø 5.5 x 18	Ø 11.5 x 2.4	Ø 22 x 3.5
Shank	104-HP	104-HP	104-HP
○ optimal / maximum rpm	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

BOBCAT F	BOBCAT126	BOBCAT164	BOBCAT154
Grit fine	Single \$53.95	Single \$28.95	Single \$34.95
			
Sizes [mm]	Ø 26 x 2	Ø 17 x 2.6	Ø 22 x 3.2
Shank	104-HP	104-HP	104-HP
○ optimal / maximum rpm	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

BOBCAT F	BOBCAT104	BOBCAT132	BOBCAT178
Grit fine	10 pack \$149.95 Single \$16.95	10 pack \$99.95 Single \$13.95	Single \$36.95
			
Sizes [mm]	Ø 5.5 x 18	Ø 11.5 x 2.4	Ø 22 x 3.5
Shank	104-HP	104-HP	104-HP
○ optimal / maximum rpm	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

LAB SUPPLIES (ROTARY)

JAGUAR

For Ceramic



Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
		■					

- Our Economical Line**
- Ideal for transition of metal and ceramic
- Economical series without diamond
- Silicon carbide
- Soft and comfortable bonding material



JAGUAR | M

Grit medium

| Sizes [mm]

| Shank

○ optimal / maximum rpm

| Mandrels

Part #	100 pack \$70.95	100 pack \$70.95	100 pack \$70.95	100 pack \$70.95
JAGUAR118	JAGUAR102	JAGUAR128	JAGUAR124	
10 pack \$8.95	10 pack \$9.95	10 pack \$9.95	10 pack \$9.95	
				
Ø 22 x 3.5	Ø 22 x 3.2	Ø 6 x 25	Ø 7 x 20	
900	900	900	900	
10,000 / 20,000	10,000 / 20,000	10,000 / 20,000	10,000 / 20,000	
A / B	A / B	F	F	

| matching mandrels please see page 54

JAGUAR | F

Grit fine

| Sizes [mm]

| Shank

○ optimal / maximum rpm

| Mandrels

Part #	100 pack \$70.95	100 pack \$70.95	100 pack \$70.95	100 pack \$70.95
JAGUAR116	JAGUAR100	JAGUAR126	JAGUAR122	
10 pack \$9.95	10 pack \$9.95	10 pack \$9.95	10 pack \$9.95	
				
Ø 22 x 3.5	Ø 22 x 3.2	Ø 6 x 25	Ø 7 x 20	
900	900	900	900	
10,000 / 20,000	10,000 / 20,000	10,000 / 20,000	10,000 / 20,000	
A / B	A / B	F	F	

| matching mandrels please see page 56

JAGUAR | XF

Grit extra fine

| Sizes [mm]

| Shank

○ optimal / maximum rpm

| Mandrels

Part #	100 pack \$70.95
JAGUAR120	JAGUAR104
10 pack \$9.95	10 pack \$9.95
	
Ø 22 x 3.5	Ø 22 x 3.2
900	900
10,000 / 20,000	10,000 / 20,000
A / B	A / B

| matching mandrels please see page 56


Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium	PEEK
								■

LEOPARD

For Peek



- ☑ **Two-step system for finishing PEEK**
- ☑ Finishing of sprue attachments on millings
- ☑ High-gloss polishing of CAD/CAM restorations
- ☑ Smooth surfaces and a perfect high-gloss luster on all pressed workpieces



LEOPARD | F

Grit fine

| Sizes [mm]


| Shank

○ optimal / maximum rpm

Part #	LEOPARD100	LEOPARD108	LEOPARD104
	10 pack \$31.99	10 pack \$53.95	10 pack \$53.95
			
	Ø 5 x 16	Ø 22 x 4	Ø 22 x 3.2
	104-HP	104-HP	104-HP
	8,000 / 20,000	8,000 / 20,000	8,000 / 20,000

LEOPARD | XF

Grit extra fine




| Sizes [mm]

| Shank

○ optimal / maximum rpm

Part #	LEOPARD102	LEOPARD110	LEOPARD106
	10 pack \$31.99	10 pack \$53.95	10 pack \$53.95
			
	Ø 5 x 16	Ø 22 x 4	Ø 22 x 3.2
	104-HP	104-HP	104-HP
	8,000 / 20,000	8,000 / 20,000	8,000 / 20,000



LAB SUPPLIES (ROTARY)

DESERT CAT

For Composite




Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
			■				

- ☑ Specially developed for veneering composites
 - ☑ New extra thin disc
- ☑ 2-step polishing system with diamond
 - ☑ Extremely fast polishing results



DESERT CAT | M




Grit medium



Sizes [mm]


| Shank

○ optimal / maximum rpm

	Single \$15.95	Single \$15.95	Single \$17.95
Part #	DESERTCAT102	DESERTCAT106	DESERTCAT110
	10 pack \$139.95	10 pack \$139.95	10 pack \$159.95
			
	Ø 5.5 x 16	Ø 15	Ø 15 x 2.3
	104-HP	104-HP	104-HP
	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

DESERT CAT | F




Grit fine



Sizes [mm]

| Shank

○ optimal / maximum rpm

	Single \$15.95	Single \$15.95	Single \$17.95
Part #	DESERTCAT100	DESERTCAT104	DESERTCAT108
	10 pack \$139.95	10 pack \$139.95	10 pack \$159.95
			
	Ø 5.5 x 16	Ø 15	Ø 15 x 2.3
	104-HP	104-HP	104-HP
	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
		○	○	○	○	■	○

LYNX
For Universal



- Our Economical Line without Diamond**
 - Smoothing with slight shine
 - Smooth binding material
 - Over 1 Million units sold



LYNX | M

Grit medium

| Sizes [mm]

| Shank 900

○ optimal / maximum rpm

| Mandrels A / B

Part #	LYNX102	LYNX104	LYNX106	LYNX108
100 pack	\$69.95	\$69.95	\$69.95	\$69.95
10 pack	\$10.95	\$10.95	\$10.95	\$10.95
Image				
Size	Ø 22 x 3.2	Ø 17 x 2.7	Ø 22 x 3.5	Ø 15 x 2.5
Shank	900	900	900	900
Optimal / max rpm	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000
Mandrels	A / B	A / B	A / B	A / B

| matching mandrels please see page 56

Part #	LYNX110	LYNX112
100 pack	\$70.95	\$70.95
10 pack	\$9.95	\$9.95
Image		
Size	Ø 7 x 20	Ø 6 x 25
Shank	900	900
Optimal / max rpm	6,000 / 20,000	6,000 / 20,000
Mandrels	F	F

| matching mandrels please see page 56



LAB SUPPLIES (ROTARY)

COUGAR

For Universal



Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
		○	○	○		■	

- ☑ **Our Economical Line with Silicium-Carbide**
 - ☑ Universal use on all dental materials
- ☑ Designed for aggressive removal of dental materials
 - ☑ Over 1 Million units sold globally



COUGAR | C

Grit coarse

Sizes [mm]

| Shank

○ optimal / maximum rpm

| Mandrels

| Sizes [mm]

| Shank










○ optimal / maximum rpm

| Mandrels

| Sizes [mm]

| Shank

○ optimal / maximum rpm

	100 pack \$70.95	100 pack \$64.95	100 pack \$64.95
Part #	COUGAR124	COUGAR132	COUGAR142
	10 pack \$9.95	10 pack \$9.95	10 pack \$9.95
			
	Ø 22 x 3.2	Ø 17 x 2.7	Ø 22 x 1
	900	900	900
	8,000 / 20,000	8,000 / 20,000	6,000 / 20,000
	A / B	A / B	A / B
	matching mandrels please see page 54		
	100 pack \$64.95	100 pack \$64.95	100 pack \$64.95
Part #	COUGAR150	COUGAR158	COUGAR174
	10 pack \$9.95	10 pack \$9.95	10 pack \$9.95
			
	Ø 22 x 4	Ø 18 x 3	Ø 6 x 25
	900	900	900
	8,000 / 20,000	8,000 / 20,000	6,000 / 20,000
	A / B	A / B	F
	matching mandrels please see page 54		
	100 pack \$20.95	100 pack \$20.95	100 pack \$20.95
Part #	COUGAR100	COUGAR108	COUGAR116
	10 pack \$20.95	10 pack \$20.95	10 pack \$20.95
			
	Ø 5 x 16	Ø 15	Ø 15 x 2.5
	104-HP	104-HP	104-HP
	6,000 / 20,000	8,000 / 20,000	8,000 / 20,000

Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
				■			

OCELOT

For Denture Acrylic



- High-grit concentration
- For removal and smoothing
- Wear-resistant

OCELOT | C

Grit: coarse

| Sizes [mm]

| Shank

○ optimal / maximum rpm

OCELOT | M

Grit: medium

| Sizes [mm]

| Shank

○ optimal / maximum rpm

Part #	OCELOT100	OCELOT104	OCELOT108
	10 pack \$37.95	10 pack \$37.95	10 pack \$41.95
			
	Ø 10 x 20.5	Ø 10 x 24	Ø 15 x 18
	104-HP	104-HP	104-HP
	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

Part #	OCELOT102	OCELOT106	OCELOT110
	10 pack \$37.95	10 pack \$37.95	10 pack \$41.95
			
	Ø 10 x 20.5	Ø 10 x 24	Ø 15 x 18
	104-HP	104-HP	104-HP
	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000



Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
				■	○	○	○

- Impregnated abrasive pad
- Extremely soft
- Efficiently applicable on many different surfaces
- Generates coarse to fine grinding pattern

PUMA | M XF UF




Grit: medium, extra fine, ultra fine

| Grit

| Sizes [mm]

| Shank

○ optimal / maximum rpm

Part #	PUMA100	PUMA104	PUMA102
	10 pack \$29.95	10 pack \$29.95	10 pack \$29.95
			
	M	XF	UF
	Ø 26	Ø 26	Ø 26
	104-HP	104-HP	104-HP
	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

PUMA

For Denture Acrylic



LAB SUPPLIES (ROTARY)

CHEETAH

For Denture Acrylic



- From removal to high-gloss
- Over 1 Million units sold
- Smooth silicone binding material

Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
				■			

CHEETAH | C

Grit coarse

| Sizes [mm]

| Shank

○ optimal / maximum rpm

Part #	CHEETAH112	CHEETAH100	CHEETAH118	CHEETAH130	CHEETAH154	CHEETAH142
	10 pack \$26.95	10 pack \$31.95	10 pack \$31.95	10 pack \$31.95	10 pack \$31.95	10 pack \$31.95
	Ø 5.5 x 16	Ø 10 x 20.5	Ø 12 x 22.5	10 x 24	Ø 7.8 x 20.5	Ø 15 x 18
	104-HP	104-HP	104-HP	104-HP	104-HP	104-HP
	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

CHEETAH | M

Grit medium

| Sizes [mm]

| Shank

○ optimal / maximum rpm

Part #	CHEETAH116	CHEETAH104	CHEETAH122	CHEETAH134	CHEETAH158	CHEETAH146
	10 pack \$26.95	10 pack \$31.95	10 pack \$31.95	10 pack \$31.95	10 pack \$31.95	10 pack \$31.95
	Ø 5.5 x 16	Ø 10 x 20.5	Ø 12 x 22.5	Ø 10 x 24	Ø 7.8 x 20.5	Ø 15 x 18
	104-HP	104-HP	104-HP	104-HP	104-HP	104-HP
	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

CHEETAH | F

Grit fine

| Sizes [mm]

| Shank

○ optimal / maximum rpm

Part #	CHEETAH114	CHEETAH102	CHEETAH120	CHEETAH132	CHEETAH156	CHEETAH144
	10 pack \$26.95	10 pack \$31.95	10 pack \$31.95	10 pack \$31.95	10 pack \$31.95	10 pack \$31.95
	Ø 5.5 x 16	Ø 10 x 20.5	Ø 12 x 22.5	Ø 10 x 24	Ø 7.8 x 20.5	Ø 15 x 18
	104-HP	104-HP	104-HP	104-HP	104-HP	104-HP
	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
					■		○

TIGER

For Metals



- Removal and smoothing of non-precious alloys**
- Wear-resistant hard binding material



TIGER | C

Grit coarse

| Sizes [mm]

| Shank

○ optimal / maximum rpm

| mandrels

Part #	TIGER118	TIGER106	TIGER114	TIGER110
	100 pack \$64.95 10 pack \$9.95	100 pack \$64.95 10 pack \$9.95	100 pack \$64.95 10 pack \$9.95	100 pack \$64.95 10 pack \$9.95
	Ø 7 x 20	Ø 22 x 3.2	Ø 22 x 4	Ø 26 x 1.3
	900	900	900	900
	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000
	F	A / B	A / B	A / B

| matching mandrels please see page 22

TIGER | M

Grit medium

| Sizes [mm]

| Shank

○ optimal / maximum rpm

| mandrels

Part #	TIGER122	TIGER124	TIGER120	TIGER108	TIGER116
	100 pack \$64.95 10 pack \$9.95	100 pack \$64.95 10 pack \$9.95	100 pack \$64.95 10 pack \$9.95	100 pack \$64.95 10 pack \$9.95	100 pack \$64.95 10 pack \$9.95
	Ø 2 x 20	Ø 3 x 24.5	Ø 7 x 20	Ø 22 x 3.2	Ø 22 x 4
	000	000	900	900	900
	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000
	G	H	F	A / B	A / B

| matching mandrels please see page 56

LAB SUPPLIES (ROTARY)

LION For Precious Metals



- Wear-resistant bonding material
- Silicium carbide
- From removal to mirror finish

	Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
LION C + M						○	■	○

LION C	Part #LION138	LION146	LION154	LION162	LION174	LION182
Grit	coarse					
	100 pack \$70.95 10 pack \$9.95	100 pack \$70.95 10 pack \$9.95	100 pack \$70.95 10 pack \$9.95	100 pack \$70.95 10 pack \$9.95	100 pack \$70.95 10 pack \$9.95	100 pack \$70.95 10 pack \$9.95
Sizes [mm]	Ø 22 x 3.2	Ø 22 x 1.2	Ø 22 x 4	Ø 7 x 20	Ø 2 x 20	Ø 3 x 24.5
Shank	900	900	900	900	000	000
○ optimal / maximum rpm	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000
Mandrels	A / B	A / B	A / B	F	G	H

| matching mandrels please see page 56

	LION98	LION106	LION114	LION122
	100 pack \$20.95	100 pack \$20.95	100 pack \$23.95	100 pack \$23.95
Sizes [mm]	Ø 5.5 x 18	Ø 5 x 16	Ø 15	Ø 15 x 2.3
Shank	104-HP	104-HP	104-HP	104-HP
○ optimal / maximum rpm	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

LION M	Part #LION142	LION150	LION158	LION166	LION178	LION186
Grit	medium					
	100 pack \$70.95 10 pack \$9.95	100 pack \$70.95 10 pack \$9.95	100 pack \$70.95 10 pack \$9.95	100 pack \$70.95 10 pack \$9.95	100 pack \$70.95 10 pack \$9.95	100 pack \$70.95 10 pack \$9.95
Sizes [mm]	Ø 22 x 3.2	Ø 22 x 1.2	Ø 22 x 4	Ø 7 x 20	Ø 2 x 20	Ø 3 x 24.5
Shank	900	900	900	900	000	000
○ optimal / maximum rpm	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000
Mandrels	A / B	A / B	A / B	F	G	H

| matching mandrels please see page 56

	LION102	LION110	LION118	LION126
	100 pack \$20.95	100 pack \$20.95	100 pack \$23.95	100 pack \$23.95
Sizes [mm]	Ø 5.5 x 18	Ø 5 x 16	Ø 15	Ø 15 x 2.3
Shank	104-HP	104-HP	104-HP	104-HP
○ optimal / maximum rpm	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

LION
For Precious Metals



Zirconia sintered	Zirconia un sintered	Ceramic	Composite	Denture acrylic	Non-precious alloys	Precious metals / Amalgam	Titanium
					■	■	

LION F + XF

LION F		100 pack \$70.95 Part #LION140 10 pack \$9.95	100 pack \$70.95 LION148 10 pack \$9.95	100 pack \$70.95 LION156 10 pack \$9.95	100 pack \$70.95 LION164 10 pack \$9.95	100 pack \$70.95 LION176 10 pack \$9.95	100 pack \$70.95 LION184 10 pack \$9.95
Grit	fine						
Sizes [mm]		Ø 22 x 3.2	Ø 22 x 1.2	Ø 22 x 4	Ø 7 x 20	Ø 2 x 20	Ø 3 x 24.5
Shank		900	900	900	900	000	000
○ optimal / maximum rpm		6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000
Mandrels		A / B	A / B	A / B	F	G	H

| matching mandrels please see page 56

	LION100 10 pack \$20.95	LION108 10 pack \$20.95	LION116 10 pack \$23.95	LION124 10 pack \$23.95
Sizes [mm]	Ø 5.5 x 18	Ø 5 x 16	Ø 15	Ø 15 x 2.3
Schaft Shank	104-HP	104-HP	104-HP	104-HP
○ optimal / maximum rpm	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

LION XF		100 pack \$70.95 Part #LION144 10 pack \$9.95	100 pack \$70.95 LION152 10 pack \$9.95	100 pack \$70.95 LION160 10 pack \$9.95	100 pack \$70.95 LION168 10 pack \$9.95	100 pack \$70.95 LION180 10 pack \$9.95	100 pack \$70.95 LION188 10 pack \$9.95
Grit	extra fine						
Sizes [mm]		Ø 22 x 3.2	Ø 22 x 1.2	Ø 22 x 4	Ø 7 x 20	Ø 2 x 20	Ø 3 x 24.5
Schaft Shank		900	900	900	900	000	000
○ optimal / maximum rpm		6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000
Mandrels		A / B	A / B	A / B	F	G	H

| matching mandrels please see page 56

	LION104 10 pack \$20.95	LION112 10 pack \$20.95	LION120 10 pack \$23.95	LION128 10 pack \$23.95
Sizes [mm]	Ø 5.5 x 18	Ø 5 x 16	Ø 15	Ø 15 x 2.3
Shank	104-HP	104-HP	104-HP	104-HP
○ optimal / maximum rpm	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000	6,000 / 20,000

PAMPAS CAT



Mandrels / Pin Holder

Mandrels and pin holders made of stainless steel for the shown items.

PAMPAS CAT | S

mandrels/pin holder

| Sizes [mm]







| Shank

○ maximum rpm

| Sizes [mm]

| Shank

○ maximum rpm

	Single \$4.95	Single \$4.95	Single \$4.95	Single \$4.95			
Part #	PAMPASCAT102	PAMPASCAT114	PAMPASCAT110	PAMPASCAT112			
	10 pack \$31.95	10 pack \$31.95	10 pack \$31.95	10 pack \$31.95			
[A]		[B]		[C]		[D]	
	Ø 5		Ø 5		Ø 8		Ø 5
	104-HP		104-HP		104-HP		104-HP
	20,000		20,000		20,000		20,000
	Single \$10.95	Single \$10.95					
	PAMPASCAT104	PAMPASCAT106					
	10 pack \$75.95	10 pack \$75.95					
[G]		[H]					
	Ø 2		Ø 3				
	104-HP		104-HP				
	20,000		20,000				




Your benefits / Vos avantages

Cutting edge technology with exceptional quality / *Qualité exceptionnelle*

- Made in Germany
Fabriqué en Allemagne
- Cutting-edge manufacturing technology
Utilisant la dernière technologie de production
- Color-coded for optimum efficiency
For best results use instruments in sequence following the colour-coding from dark to light (1-6).
Carbides and diamonds are interchangeable as long as you follow the color-coding.

Code couleur pour une efficacité optimale

Pour obtenir rapidement d'excellents résultats, utilisez les instruments dans l'ordre du code couleur. Les fraises carbures et diamantées sont interchangeables à condition que vous suiviez le code couleur du plus foncé au plus clair (1-6).



- Note: We recommend not to skip more than one color.
Moving from black to blue (1-3) is possible; moving from black directly to yellow (1-5) is inefficient.
Remarque : Nous vous conseillons de ne pas sauter plus d'une couleur.
Exemple : 1 à 3 est possible, mais 1 à 5 n'est pas recommandé, car inefficace et/ou perte de temps.
- Each instrument is color-coded according to cut/grit size.
Les bagues couleur correspondent à la denture / au grain et se trouvent sur la tige de l'instrument.

DIADUR the new generation of carbide cutters

DIADUR: la nouvelle génération de fraises carbure de tungstène


The new micro - structure carbide cutter is manufactured with today's most advanced technology.

Les nouvelles fraises carbure de tungstène microstructures sont fabriquées à l'aide d'une technologie de pointe

Benefits and advantages


Vos avantages

- Extended lifetime even in case of intensive use
Durée de vie prolongée même en cas d'utilisation intensive
- Precise instrument balance (less than 0.015 mm), therefore smoother running of your handpiece, reduces vibration and protects chuck and bearings
La rotation d'une précision exceptionnellement élevée (< 0,015 mm) permet un fonctionnement silencieux sans vibrations tout en protégeant vos pièces à main
- Optimal surface of workpieces
Qualité de surface optimale de vos pièces de travail
- Improved cutting efficiency due to innovative cutting geometries
Abrasion améliorée grâce à des géométries de coupe innovantes
- Easy to work even after extended use
Travail facile même après une utilisation de longue durée



Unlike other manufacturers, we solder the shanks into the head of our instruments resulting in improved balance and safety.


Pour des raisons de sécurité et pour améliorer la précision de rotation, la tige est brasée dans la tête de tungstène.



Uni-disc

Uni-disc

Keeps chuck and bearings clean
Maintien propre de la pince de serrage et des roulements



Correct insertion of instrument!
Insertion correcte de l'instrument !



Always wear safety goggles!
Toujours porter des lunettes de protection !

LAB SUPPLIES (ROTARY)

DIADUR carbide cutters / DIADUR fraises en carbure de tungstène

TURBO



15' - 25.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

Price in CAD



extracoarse, black ring / extragros, bague noire

060 1	040 1	060 1	070 1	060 1	060 1
509901	509201	557701	509101	509001	559001
●	●	●	●	●	●
22.99	20.99	22.99	22.99	17.99	23.99

Application/Utilisation



- Plaster
Plâtre
- Thermoplastics
Thermoplastiques
- Tray material
Matériel d'empreinte
- Epoxy models
Modèles époxy
- Heat and self curing acrylics
Acryliques pour prothèses

TRI-CUTTER + CROSS-CUTTER



5' - 10.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

Price in CAD



extracoarse, black ring / extragros, bague noire

023 1	023 1
80049	80055
●	●
13.99	13.99



- TRI-CUTTER
- For cutting of thermoplastic and tray material.
Pour couper tous les matériaux thermoplastiques et d'empreinte
- CROSS-CUTTER
- Acrylics and plaster
Acryliques et plâtre

MACRO




15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

Price in CAD



coarse, green ring / gros, bague verte

023 1	040 1	060 1	023 1	023 1	023 1	060 1
406702	405001	408301	408502	406502	406402	407501
●	●	●	●	●	●	●
17.99	17.99	17.99	17.99	18.99	17.99	20.99



- Cobalt chrome alloys
Alliages Co-Cr
- Plaster
Plâtre
- Heat, self and light curing acrylics
Acryliques pour prothèses

DIADUR carbide cutters / DIADUR fraises en carbure de tungstène

MACRO



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE



coarse, green ring / gros, bague verte

023	023	040	045	050	060	070
1	1	1	1	1	1	1

A06302	A06102	A06801	A03201	A06901	A07201	A07401
--------	--------	--------	--------	--------	--------	--------

Order-Nr.

Price in CAD

17.99	18.99	20.99	21.99	21.99	23.99	19.99
-------	-------	-------	-------	-------	-------	-------

Application/Utilisation



- Cobalt chrome alloys
Alliages Co-Cr
- Plaster
Plâtre
- Heat, self and light curing acrylics
Acryliques pour prothèses

MACRO



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE



coarse, green ring / gros, bague verte

060	070	040	060	060	023	060
1	1	1	1	1	1	1

A08101	A07101	A03301	A07001	A07301	A06602	A07601
--------	--------	--------	--------	--------	--------	--------

Order-Nr.

Price in CAD

17.99	26.99	21.99	17.99	27.99	17.99	24.99
-------	-------	-------	-------	-------	-------	-------



- Cobalt chrome alloys
Alliages Co-Cr
- Plaster
Plâtre
- Heat, self and light curing acrylics
Acryliques pour prothèses

QUATTRO



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE



medium, blue ring / standard, bague bleue

023	040	050	060	023	023	023	016
1	1	1	1	1	1	1	1

606702	605001	605201	607901	601402	602102	602302	605902
--------	--------	--------	--------	--------	--------	--------	--------

Order-Nr.

Price in CAD

17.99	21.99	23.99	28.99	17.99	17.99	17.99	18.99
-------	-------	-------	-------	-------	-------	-------	-------



- Cobalt chrome alloys
Alliages Co-Cr
- Non-precious alloys
Alliages non-précieus
- Precious alloys
Alliages précieux
- Semi-precious alloys
Alliages semi-précieus

DIADUR carbide cutters / DIADUR fraises en carbure de tungstène

QUATTRO



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE



medium, blue ring / standard, bague bleue

023	031	040	045	060	050	060	060
1	1	1	1	1	1	1	1

606102	605502	606801	603201	602701	602901	607201	607501
--------	--------	--------	--------	--------	--------	--------	--------

Order-Nr.

Price in CAD

17.99	19.99	20.99	21.99	20.99	24.99	25.99	25.99
-------	-------	-------	-------	-------	-------	-------	-------

Application/Utilisation



- Cobalt chrome alloys
Alliages Co-Cr
- Non-precious alloys
Alliages non-précieux
- Precious alloys
Alliages précieux
- Semi-precious alloys
Alliages semi-précieux

QUATTRO



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE



medium, blue ring / standard, bague bleue

023	029	060	060	040	060	023	060
1	1	1	1	1	1	1	1

604702	605302	608101	603401	603301	602001	606602	607601
--------	--------	--------	--------	--------	--------	--------	--------

Order-Nr.

Price in CAD

16.99	19.99	28.99	30.99	20.99	18.99	18.99	26.99
-------	-------	-------	-------	-------	-------	-------	-------



- Cobalt chrome alloys
Alliages Co-Cr
- Non-precious alloys
Alliages non-précieux
- Precious alloys
Alliages précieux
- Semi-precious alloys
Alliages semi-précieux

STRATO



15' - 25.000 rpm/tpm

ø 1/10 mm

Pu / UE



medium, blue ring / standard, bague bleue

023	040	040	040	060
1	1	1	1	1

416102	616801	415401	413301	417001
--------	--------	--------	--------	--------

Order-Nr.

Price in CAD

18.99	20.99	20.99	20.99	22.99
-------	-------	-------	-------	-------



- Heat, self and light curing acrylics
Acryliques pour prothèses
- Soft lining material
Acryliques mous
- Valplast and similar
Valplast et similaires

DIADUR carbide cutters / DIADUR fraises en carbure de tungstène

QUATTRO-TWIN




15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

Price in CAD



medium, blue ring / standard, bague bleue

023
1

602103



28.99

Application/Utilisation



- Cobalt chrome alloys
Alliages Co-Cr
- Non-precious alloys
Alliages non-précieus
- Precious alloys
Alliages précieux
- Semi-precious alloys
Alliages semi-précieus
- no uni-disc
sans Uni-disc

MACRONOM



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

Price in CAD



coarse, green ring / gros, bague verte

060
1

422001



20.99



- Plaster
Plâtre
- Acrylics
Acryliques

ERGONOM



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

Price in CAD



fine, red ring / fin, bague rouge

023
1

023
1

040
1

060
1

100302

100202

100801

101001



17.99

17.99

19.99

21.99



- Precious alloys
Alliages précieux

DIADUR carbide cutters / DIADUR fraises en carbure de tungstène

MICRO



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Order-Nr.

Price in CAD

fine, red ring / fin, bague rouge							
023	023	060	023	023	045	018	023
1	1	1	1	1	1	1	1
302502	305102	303901	302402	302302	302802	304902	304802
●	●	●	●	●	●	●	●
19.99	19.99	28.99	19.99	19.99	22.99	20.99	20.99

Application/Utilisation



- Precious alloys
Alliages précieux
- Semi-precious alloys
Alliages semi-précieux
- Cobalt chrome alloys
Alliages Co-Cr
- Non-precious alloys
Alliages non-précieux
- Composites
Composites

MICRO



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Order-Nr.

Price in CAD

fine, red ring / fin, bague rouge							
023	023	031	040	045	050	023	023
1	1	1	1	1	1	1	1
302202	306102	305502	302801	303201	302901	304702	304602
●	●	●	●	●	●	●	●
19.99	21.99	23.99	22.99	23.99	24.99	19.99	19.99



- Precious alloys
Alliages précieux
- Semi-precious alloys
Alliages semi-précieux
- Cobalt chrome alloys
Alliages Co-Cr
- Non-precious alloys
Alliages non-précieux
- Composites
Composites

MICRO



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Order-Nr.

Price in CAD

fine, red ring / fin, bague rouge						
060	040	060	060	023	023	060
1	1	1	1	1	1	1
308101	303301	303001	307301	305002	302602	303601
●	●	●	●	●	●	●
27.99	22.99	25.99	27.99	18.99	18.99	25.99



- Precious alloys
Alliages précieux
- Semi-precious alloys
Alliages semi-précieux
- Cobalt chrome alloys
Alliages Co-Cr
- Non-precious alloys
Alliages non-précieux
- Composites
Composites

DIADUR carbide cutters / DIADUR fraises en carbure de tungstène

MATRO



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

Price in CAD


fine, red ring / fin, bague rouge						
023	060	016	023	040	040	060
1	1	1	1	1	1	1
315102	313901	314902	314802	312801	313301	313001
●	●	●	●	●	●	●
17.99	19.99	17.99	17.99	18.99	18.99	19.99

Application/Utilisation



- Heat, self and light curing acrylics
Acryliques pour prothèses
- Soft lining material
Acryliques mous
- Valplast and similar
Valplast et similaires

MILLI-MICRON



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

Price in CAD

superfine, yellow ring / superfin, bague jaune							
023	016	023	045	023	060	023	023
1	1	1	1	1	1	1	1
201402	204902	201202	203201	204602	203001	205002	201602
●	●	●	●	●	●	●	●
20.99	21.99	19.99	26.99	21.99	27.99	23.99	19.99



- Non-precious alloys
Alliages non-précieux
- Semi-precious alloys
Alliages précieux
- Precious alloys
Alliages semi-précieux
- Composites
Composites

SUPER-TECH



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

Price in CAD

fine, red ring / fin, bague rouge	
040	1
902801	●
30.99	



- Titanium
Titane
- Cobalt chrome alloys
Alliages Co-Cr
- Non-precious alloys
Alliages non-précieux

**3-way sectional cut
Triple denture**

DIADUR carbide cutters / DIADUR fraises en carbure de tungstène

TITANO



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Order-Nr.

Price in CAD

extracoarse, black ring / extragros, bague noire



023
1

736602F



23.99

Application/Utilisation



- Titanium
Titane

TITANO



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Order-Nr.

Price in CAD

fine, red ring / fin, bague rouge



023
1

702402F



22.99



023
1

702302F



23.99



010
1

101102F




9.99



- Titanium
Titane
- Non-precious alloys
Alliages non-précieux
- Composites
Composites
- Light curing tray material
Matériaux de prise d'empreinte
photopolymérisables

TITANO



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Order-Nr.

Price in CAD

superfine, yellow ring / superfin, bague jaune



040
1

7008015F



29.99



- Titanium
Titane
- Non-precious alloys
Alliages non-précieux
- Composites
Composites
- Light curing tray material
Matériaux de prise d'empreinte
photopolymérisables

DIADUR carbide cutters / DIADUR fraises en carbure de tungstène

STAR



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

Price in CAD

	ultrafine, white ring / ultrafin, bague blanche	
	040	060
	1	1
	700801	708001
	○	○
	22.99	25.99

Application/Utilisation



- Precious alloys
Alliages précieux
- Tough acrylic veneers
Acryliques facettes extra durs
- Titanium
Titane

SUPER-STAR



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

Price in CAD

	ultrafine, white ring / ultrafin, bague blanche	
	023	
	1	
	712602SF	
	○	
	23.99	



- Non-precious alloys
Alliages non-précieux
- Composites
Composites
- Soft lining material
Acryliques mous

PINBUR



ø 1/10 mm
Pu / UE

Order-Nr.

Price in CAD

Blade length <i>Longueur des lames</i>	5,5 mm { 13,3 mm {	{ d ₁ { { d ₂ {	
	without ring / sans bague		
	d ₁ = 2mm / d ₂ = 3mm		
	1		
	80069		
	●		
	27.99		

- Bur for PINDEX® system
Fraise pour système PINDEX®
- Total tool length 35 mm
Longueur totale de l'instrument 35 mm

PINDEX® is a registered trademark of Coltene/Whaledent, Inc USA
PINDEX® est une marque déposée de Coltene/Whaledent, Inc USA

DIADUR carbide instruments / *DIADUR instruments en carbure de tungstène*

1 Rosehead burs / *Fraises rondes*



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

US-Code

Price in CAD

005 6	006 6	008 6	010 6	012 6	014 6	016 6	018 6	021 6	023 6	027 6
117005	117006	117008	117010	117012	117014	117016	117018	117021	117023	117027
●	●	●	●	●	●	●	●	●	●	●
1/4	1/2	1	2	3	4	5	6	7	8	10
34.99	34.99	34.99	34.99	34.99	34.99	34.99	34.99	34.99	34.99	50.99

2 Inverted cone burs / *Fraises coniques renversées*



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

US-Code

Price in CAD

008 6	010 6	014 6	016 6
113008	113010	113014	113016
●	●	●	●
34	35	37	38
34.99	34.99	34.99	34.99

23



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

US-Code

Price in CAD

008 6	010 6
108608	108610
●	●
168	170
34.99	34.99

23L



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

US-Code

Price in CAD

009 6
108709
●
169L
34.99

23R



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.


US-Code

Price in CAD





008 6	010 6
101101	101102
●	●
1169	1170
34.99	34.99

DIADUR carbide instruments / DIADUR instruments en carbure de tungstène

31




15' - 30.000 rpm/tpm
ø 1/10 mm
Pu / UE

			
010 6	014 6	016 6	023 6
118010	118014	118016	118023
●	●	●	●
557 34.99	559 34.99	560 34.99	563 34.99



Order-Nr.

US-Code
Price in CAD

31L




15' - 30.000 rpm/tpm
ø 1/10 mm
Pu / UE

	
010 6	012 6
108810	108812
●	●
557L 34.99	558L 34.99





Order-Nr.

US-Code
Price in CAD

33




15' - 30.000 rpm/tpm
ø 1/10 mm
Pu / UE

			
008 6	012 6	016 6	021 6
109308	109312	109316	109321
●	●	●	●
699 34.99	701 34.99	702 34.99	703 34.99



Order-Nr.

US-Code
Price in CAD

33L




15' - 30.000 rpm/tpm
ø 1/10 mm
Pu / UE

	
008 6	010 6
109408	109410
●	●
699L 34.99	700L 34.99


Order-Nr.

US-Code
Price in CAD

OKKLU-FIN




15' - 30.000 rpm/tpm
ø 1/10 mm
Pu / UE



023 6
304402
●
68.99

Order-Nr.

Price in CAD



15' - 30.000 rpm/tpm
ø 1/10 mm
Pu / UE



014 6
118080
●
62.99

Order-Nr.

Price in CAD

SIDIA sintered diamonds / *SIDIA Diamants frittés*

FG MACRO



Included free in each package!
Y compris réducteur !

30' - 300.000 rpm/tpm

ø 1/10 mm
Pu / UE

medium, blue ring / standard, bague bleue

016
1

56212

Order-Nr.

Price in CAD

23.99

Application/Utilisation



- Ceramics
Céramique

HP MACRO

15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

medium, blue ring / standard, bague bleue

016
1

58242

Order-Nr.

Price in CAD

23.99

016
1

58202

23.99



- Ceramics
Céramique

SUPERTURBO

15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

extra coarse, black ring / extragros, bague noire

050
1

53021

Order-Nr.

Price in CAD

39.99

060
1

53161

39.99



- Cobalt chrome alloys
Alliages Co-Cr

SIDIA sintered diamonds / SIDIA Diamants frittés

TURBO

15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

	coarse, green ring / gros, bague verte		
	050	050	060
	1	1	1
	55021	55181	55161
Order-Nr.	●	●	●
Price in CAD	39.99	34.99	39.99

Application/Utilisation



- Cobalt chrome alloys
Alliages Co-Cr
- Non-precious alloys
Alliages non-précieus
- Semi-precious alloys
Alliages semi-précieus

MACRO

15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

	medium, blue ring / standard, bague bleue
	060
	1
	59202
Order-Nr.	●
Price in CAD	29.99



- Cobalt chrome alloys
Alliages Co-Cr
- Ceramics
Céramique

MICRO

15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

	fine, red ring / fin, bague rouge						
	060	080	027	027	030	025	025
	1	1	1	1	1	1	1
	50650	50051	50215	50192	50082	50212	50217
Order-Nr.	●	●	●	●	●	●	●
Price in CAD	37.99	34.99	32.99	29.99	32.99	32.99	29.99



- Cobalt chrome alloys
Alliages Co-Cr
- Ceramics
Céramique

LAB SUPPLIES (ROTARY)

SIDIA sintered diamonds / *SIDIA Diamants frittés*

MICRO

15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Order-Nr.

Price in CAD


	fine, red ring / fin, bague rouge						
	023	040	030	030	030	060	080
	1	1	1	1	1	1	1
	50202	50044	50172	50092	50214	50161	50131
	●	●	●	●	●	●	●
	32.99	34.99	34.99	32.99	32.99	39.99	32.99

Application/Utilisation



- Cobalt chrome alloys
Alliages Co-Cr
- Ceramics
Céramique

CLEANING STONE / *PIERRE DE NETTOYAGE*




Pu / UE

Order-Nr.

Price in CAD


	1
	15001
	●
	7.99



- Cleaning stone for sintered and electroplated diamonds
Pour le nettoyage d'outils diamantés frittés et galvaniques
- Apply water before using
Mouiller avant usage

SIDIA sintered diamond discs / SIDIA Disques diamantés frittés

MACRO



5' - 15.000 rpm/tpm

ø 1/10 mm
Pu / UE

medium, blue ring / standard, bague bleue

220
1

58154

Order-Nr.

Price in CAD


23.99

Application/Utilisation



- For contouring of CoCr-alloys
Pour le contourage d'alliages Co-Cr
- All sintered diamond discs are mounted and centered
Tous les disques en diamanté fritté sont montés et centrés

MICRO



5' - 15.000 rpm/tpm

ø 1/10 mm
Pu / UE

fine, red ring / fin, bague rouge

160
1

50265

Order-Nr.

Price in CAD

23.99



- For contouring of CoCr-alloys
Pour le contourage d'alliages Co-Cr
- All sintered diamond discs are mounted and centered
Tous les disques en diamanté fritté sont montés et centrés

Zirconia instruments / Instruments pour zircone

MINI-TRI-ZIR



5' - 10.000 rpm/tpm

ø 1/10 mm

Pu / UE

extracoarse / extragros

016

1

80059Zir

Order-Nr.

Price in CAD

13.99

Application/Utilisation




MINI-TRI-ZIR

- Cutting of connecting bars in green state
- Coupe des barres de liaison du zircone non-fritté

No water cooling needed!
Travail à sec!

SIDIAFLEX GOLD



3' - 5.000 rpm/tpm

ø 1/10 mm

Pu / UE

superfine, yellow ring / superfin, bague jaune

320


1

181220

Order-Nr.

Price in CAD

55.99



- Cutting of connecting bars
- Séparation des barres de liaison

Water cooling!
Refroidissement à eau!

Electroplated spiral diamond finishers / Fraises diamantées avec spirale

**PERIPHERY FINISHER
ALVEOLAR FINISHER
FINISSEUR PÉRIPHÉRIQUE
FINISSEUR ALVÉOLAIRE**

5' - 15.000 rpm/tpm

ø 1/10 mm
Pu / UE

extracoarse, black ring/extragros, bague noire

060
1

060
1

60016

60015

Order-Nr.

Price in CAD

16.99

22.99

Application/Utilisation



- Acrylics
Acryliques
- Light curing tray materials
Matériaux de prise d'empreinte photopolymérisables

PINFIX

5' - 15.000 rpm/tpm

ø 1/10 mm
Pu / UE

extracoarse, black ring/extragros, bague noire


053
1

60017

Order-Nr.

Price in CAD

22.99



• Plaster
Plâtre

**TURBO DIAMOND GRINDERS
TURBO MEULETTE DIAMANTÉE**

3' - 15.000 rpm/tpm

ø 1/10 mm
Pu / UE

extracoarse, black ring/extragros, bague noire

027
1

027
1

60023

60021

Order-Nr.

Price in CAD

14.99

14.99



- Ideal for soft reline materials
Idéal pour les matériaux de rebasage mous
- All acrylics
Tous acryliques
- Plaster
Plâtre
- Light curing tray material
Matériaux de prise d'empreinte photopolymérisables

Electroplated HP diamonds / *Fraises diamantées galvaniques PM*

801



5' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Grit size / Grain

● medium / standard

Order-Nr.

Price in CAD

009 3	010 3	012 3	014 3	016 3	018 3	021 3	023 3	029 3	035 3	040 3
60001098	60001108	60001128	60001148	60001168	60001188	60001218	60001238	60001298	60001358	60001408
●	●	●	●	●	●	●	●	●	●	●
19.99	19.99	19.99	19.99	19.99	19.99	19.99	19.99	19.99	19.99	19.99

805



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Grit size / Grain

● medium / standard

Order-Nr.

Price in CAD

014 3	018 3	025 3
60010148	60010188	60010258
●	●	●
19.99	19.99	19.99

085



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Grit size / Grain

● medium / standard

Order-Nr.

Price in CAD

023 3
60014238
●
19.99

Application/Utilisation



- Ceramics
Céramique
- Hard composites
Composites durs

836



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Grit size / Grain

● coarse / gros

Order-Nr.

Price in CAD

049 3
60110495
●
19.99

837



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Grit size / Grain





● medium / standard






Order-Nr.





Price in CAD

014 3	016 3	033 3	040 3	070 3
60111148	60111168	60111338	60111408	60111708
●	●	●	●	●
19.99	19.99	19.99	19.99	19.99

Electroplated HP diamonds / Fraises diamantées galvaniques PM

<p>837L</p>  <p>15' - 30.000 rpm/tpm</p> <p>ø 1/10 mm Pu / UE</p> <p>Grit size / Grain ● medium / standard</p> <p>Order-Nr.</p> <p>Price in CAD</p>	 <p>018 3</p> <p>027 3</p> <p>050 3</p> <p>60112188</p> <p>60112278</p> <p>60112508</p> <p>19.99 19.99 19.99</p>	<p>881</p>  <p>15' - 30.000 rpm/tpm</p> <p>ø 1/10 mm Pu / UE</p> <p>Grit size / Grain ● medium / standard ● coarse / gros</p> <p>Order-Nr.</p> <p>Price in CAD</p>	 <p>016 3</p> <p>018 3</p> <p>018 3</p> <p>025 3</p> <p>60141168</p> <p>60141188</p> <p>60141185</p> <p>60141258</p> <p>19.99 19.99 19.99 19.99</p>
--	---	---	---

<p>882L</p>  <p>15' - 30.000 rpm/tpm</p> <p>ø 1/10 mm Pu / UE</p> <p>Grit size / Grain ● medium / standard ● extracoarse / extragros</p> <p>Order-Nr.</p> <p>Price in CAD</p>	 <p>018 3</p> <p>018 3</p> <p>60143188</p> <p>60143183</p> <p>19.99 19.99</p>	<p>837KR</p>  <p>15' - 30.000 rpm/tpm</p> <p>ø 1/10 mm Pu / UE</p> <p>Grit size / Grain ● medium / standard ● coarse / gros</p> <p>Order-Nr.</p> <p>Price in CAD</p>	 <p>055 3</p> <p>055 3</p> <p>60158558</p> <p>60158555</p> <p>19.99 19.99</p>	<p>Application/Utilisation</p>  <ul style="list-style-type: none"> • Ceramics Céramique • Hard composites Composites durs <p>Black, extracoarse / Noir, extragros</p> <ul style="list-style-type: none"> • Co-Cr alloys Alliages non-précieux
---	---	--	--	---

<p>858</p>  <p>15' - 30.000 rpm/tpm</p> <p>ø 1/10 mm Pu / UE</p> <p>Grit size / Grain ● medium / standard</p> <p>Order-Nr.</p> <p>Price in CAD</p>	 <p>012 3</p> <p>60165128</p> <p>19.99</p>	<p>859</p>  <p>15' - 30.000 rpm/tpm</p> <p>ø 1/10 mm Pu / UE</p> <p>Grit size / Grain ● medium / standard</p> <p>Order-Nr.</p> <p>Price in CAD</p>	 <p>012 3</p> <p>014 3</p> <p>016 3</p> <p>018 3</p> <p>60166128</p> <p>60166148</p> <p>60166168</p> <p>60166188</p> <p>19.99 19.99 19.99 19.99</p>
---	---	---	---

Electroplated HP diamonds / Fraises diamantées galvaniques PM

859L



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Grit size / Grain
medium / standard

Order-Nr.

Price in CAD

014
3

60167148



19.99

845



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Grit size / Grain
medium / standard

Order-Nr.

Price in CAD

014
3

60170148



19.99

Application/Utilisation



- Ceramics
Céramique
- Hard composites
Composites durs

Black, extra coarse / Noir, extragros
• Co-Cr alloys
Alliages non-précieux



847



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Grit size / Grain
 medium / standard
 extra coarse / extragros

Order-Nr.

Price in CAD

018
3

60172188



19.99

033
3

60172333



19.99

040
3

60172408



19.99

050
3

60172508



19.99


850



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Grit size / Grain
 medium / standard

Order-Nr.

Price in CAD

014
3

60199148



19.99

016
3

60199168



19.99

025
3

60199258



19.99


807



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Grit size / Grain
 medium / standard

Order-Nr.

Price in CAD

018
3

60225188



19.99

023
3

60225238



19.99



860



15' - 30.000 rpm/tpm

ø 1/10 mm

Pu / UE

Grit size / Grain
 fine / fin
 medium / standard

Order-Nr.

Price in CAD

012
3

60245126



19.99

014
3

60245146



19.99

014
3


60245148



19.99

Electroplated HP diamonds / Fraises diamantées galvaniques PM

862



15' - 30.000 rpm/tpm


ø 1/10 mm
Pu / UE

Grit size / Grain
● medium / standard

Order-Nr. ●


Price in CAD 19.99

60249188



018
3

863



15' - 30.000 rpm/tpm


ø 1/10 mm
Pu / UE

Grit size / Grain
● fine / fin
● medium / standard

Order-Nr. ● ● ● ●


Price in CAD 19.99 19.99 19.99 19.99

60250108 60250128 60250166 60250168



010 012 016 016
3 3 3 3

368



15' - 30.000 rpm/tpm


ø 1/10 mm
Pu / UE

Grit size / Grain
● medium / standard

Order-Nr. ●


Price in CAD 19.99

60254478



047
3

368



15' - 30.000 rpm/tpm


ø 1/10 mm
Pu / UE

Grit size / Grain
● fine / fin
● medium / standard

Order-Nr. ● ●

Price in CAD 19.99 19.99

60257236 60257238




023 023
3 3

Application/Utilisation



369



15' - 30.000 rpm/tpm


ø 1/10 mm
Pu / UE

Grit size / Grain
● medium / standard

Order-Nr. ●


Price in CAD 19.99

60263608



060
3

825



15' - 30.000 rpm/tpm


ø 1/10 mm
Pu / UE

Grit size / Grain
● medium / standard

Order-Nr. ● ●

Price in CAD 19.99 19.99

60304408 60304608




040 060
3 3



Electroplated HP diamonds / *Fraises diamantées galvaniques PM*

889



15' - 30.000 rpm/tpm

ø 1/10 mm
Pu / UE

Grit size / Grain
● medium / standard

Order-Nr.

Price in CAD

010
3

60540108

●

19.99

Application/Utilisation



- Ceramics
Céramique
- Hard composites
Composites durs

Electroplated sintered diamond disc / *Disque galvaniques frittés*

SIDIA-FLEX

Flexible disc with true sintered diamond edge!


Disque flexible avec bord diamanté fritté!

3' - 5.000 rpm/tpm

ø 1/10 mm
Pu / UE

Order-Nr.

Price in CAD



superfine, yellow ring / superfin, bague jaune

200
1

182021

●

53.99



- For separating of:
- Ceramics incl. zirconia
 - Extremely long lasting


Pour couper et séparer :

- Céramiques y compris le zircon
- Durée de vie utile exceptionnellement élevée

Flexible diamond discs / Disques diamantés flexibles

SPIROFLEX

double sided
double-face



5' - 10.000 rpm/tpm

ø 1/10 mm
Pu / UE

fine, red ring / fin, bague rouge

220
1

122291

Order-Nr.



Price in CAD

18.99


Application/Utilisation



- Ceramics
Céramique
- Composites
Composites
- Acrylics
Acryliques
- mounted and centered
monté et centré

PLEXOFLEX

double sided
double-face



3' - 5.000 rpm/tpm

ø 1/10 mm
Pu / UE

fine, red ring / fin, bague rouge

220
1

162221

Order-Nr.



Price in CAD

20.99

220
1

192221




20.99



- For delicate separation and contouring
Pour la séparation et le contourage fin
- Ceramics
Céramique
- Tough acrylic veneers
Acryliques facettes extra durs
- mounted and centered
monté et centré

SEPAFLEX

double sided
double-face



3' - 5.000 rpm/tpm

ø 1/10 mm
Pu / UE

fine, red ring / fin, bague rouge

220
1

182221

Order-Nr.



Price in CAD

20.99

220
1

132221



20.99

220
20

13222-18



249.99




- Ceramics
Céramique
- Composites
Composites
- mounted and centered
monté et centré
- 20 pcs incl. 2 mandrels
20 non-montés plus 2 mandrins

Flexible diamond discs / Disques diamantés flexibles

POLYFLEX | CIRCO-FLEX

double sided
double-face



3' - 5.000 rpm/tpm

ø 1/10 mm
Pu / UE

fine, red ring / fin, bague rouge

220
1
182231

220
1
182251

Order-Nr.


Price in CAD

19.99

19.99

TECHNO

double sided
double-face



3' - 5.000 rpm/tpm

ø 1/10 mm
Pu / UE

fine, red ring / fin, bague rouge

220
1
187221


Order-Nr.

Price in CAD

19.99

PERFOCUT

double sided
double-face



3' - 5.000 rpm/tpm

ø 1/10 mm
Pu / UE

fine, red ring / fin, bague rouge

210
1
272121

210
1
272125

Order-Nr.

Price in CAD

18.99

18.99

Application/Utilisation



- Ceramics
Céramique
- Composites
Composites
- mounted and centered
monté et centré
- 20 pcs incl. 2 mandrels
20 non-montés plus 2 mandrins



- Ceramics
Céramique
- Composites
Composites
- mounted and centered
monté et centré




- For separating
Pour séparation
- Ceramics
Céramique
- Plaster
Plâtre
- Acrylics
Acryliques
- Tough acrylic veneers
Acryliques facettes extra durs
- mounted and centered
monté et centré

Flexible diamond discs / Disques diamantés flexibles

CUTFLEX

double sided
double-face



3' - 5.000 rpm/tpm

ø 1/10 mm
Pu / UE

fine, red ring / fin, bague rouge

220
1

152211

Order-Nr.

Price in CAD

19.99


Application/Utilisation



- For separating
Pour séparation
- Ceramics
Céramique
- Tough acrylic veneers
Acryliques facettes extra durs
- mounted and centered
monté et centré

TOTO-FLEX

double sided
double-face



5' - 10.000 rpm/tpm

ø 1/10 mm
Pu / UE

fine, red ring / fin, bague rouge

220
1

112233

Order-Nr.

Price in CAD


17.99



- For separation and contouring
Pour séparer et contourer
- Ceramics
Céramique
- Plaster
Plâtre
- Acrylics
Acryliques
- mounted and centered
monté et centré

SUPER-FLEXIFLEX

double sided
double-face



5' - 10.000 rpm/tpm

ø 1/10 mm
Pu / UE

fine, red ring / fin, bague rouge

190
1

220
1

220
1

220
20

121921

122214

122221

122218

Order-Nr.

Price in CAD

16.99

17.99

17.99

169.99




- For separation and contouring
Pour séparer et contourer
- Ceramics
Céramique
- mounted and centered
monté et centré
- 20 pcs incl. 2 mandrels
20 non-montés plus 2 mandrins

Flexible diamond discs / Disques diamantés flexibles

TOTO-FLEX

double sided
double-face



5' - 10.000 rpm/tpm

ø 1/10 mm
Pu / UE

coarse, green ring / gros, bague verte

220	220	220
1	20	1
112231	1122318	112241
●	●	●
17.99	219.99	17.99

Order-Nr.

Price in CAD

Application/Utilisation



- For separation and contouring
Pour séparer et contourer
- Ceramics
Céramique
- Plaster
Plâtre
- Acrylics
Acryliques
- mounted and centered
monté et centré
- 20 pcs incl. 2 mandrels
20 non-montés plus 2 mandrins

CONO-DISC

rigid, not flexible
rigide, non flexible



5' - 15.000 rpm/tpm

ø 1/10 mm
Pu / UE

coarse, green ring / gros, bague verte

190
1
161911

Order-Nr.

Price in CAD


25.99



- Grinding and contouring of ceramics and composites
Finition et contourage de la céramique et des composites

TORO-FLEX

double sided
double-face



5' - 10.000 rpm/tpm

ø 1/10 mm
Pu / UE

coarse, green ring / gros, bague verte

220
1

112242

Order-Nr.

Price in CAD

17.99




- For separation and contouring
Pour la séparation et le meulage latéral
- Ceramics
Céramique
- mounted and centered
monté et centré

Flexible diamond discs / Disques diamantés flexibles

FLEXIFLEX

double sided
double-face



5' - 10.000 rpm/tpm

ø 1/10 mm
Pu / UE

coarse, green ring / gros, bague verte

190	190	220	220	220
1	20	1	1	20
111921	1119218	112214	112221	1122218

Order-Nr.

Price in CAD

16.99 189.99 17.99 16.99 189.99


Application/Utilisation



- For separation and contouring
Pour la séparation et le meulage latéral
- Ceramics
Céramique
- mounted and centered
monté et centré
- 20 pcs incl. 2 mandrels
20 non-montés plus 2 mandrins

TOP-FLEX | MINI-FLEX

double sided
double-face



5' - 10.000 rpm/tpm

ø 1/10 mm
Pu / UE

superfine, yellow ring / superfin, bague jaune

190	220	090
1	1	1
141921	142221	140921

Order-Nr.

Price in CAD

17.99 19.99 14.99



- For separation and contouring
Pour séparer et contourer
- Ceramics
Céramique
- mounted and centered
monté et centré

HYDROFLEX

double sided
double-face



3' - 5.000 rpm/tpm

ø 1/10 mm
Pu / UE

ultrafine, white ring / ultrafin, bague blanche

190	190
1	1
131961	131971

Order-Nr.

Price in CAD


24.99 20.99



- Only for finishing
Uniquement pour finition
- Ceramics
Céramique
- mounted and centered
monté et centré

Separating discs & plaster discs / *Disques séparateurs & disques plâtre*

SHARK DISCS



.... rpm/tpm →

ø 1/10 mm

Pu / UE

5' - 30.000 5' - 50.000

400

260

20

20

1334018

1332638

Order-Nr.

Price in CAD

110.99


77.99

Application/Utilisation



- Reinforced cutting discs for non-precious metals with exceptional long lifetime
- Disques renforcés pour tronçonner alliages non-précieux, durée de vie exceptionnelle*

PLASTER DISCS MINI & MAXI



5' - 10.000 rpm/tpm

ø 1/10 mm

Pu / UE

fine, red ring / fin, bague rouge

320

400

1

1

153222

154022


Order-Nr.

Price in CAD

46.99

52.99

HIGH-TECH PLASTRO



5' - 10.000 rpm/tpm

ø 1/10 mm

Pu / UE

fine, red ring / fin, bague rouge

350

370

450

1

1

1

153522

153722

154522

Order-Nr.

Price in CAD

59.99

59.99

69.99



- Plaster
- Plâtre*
- Synthetic resins
- Acryliques*

Plaster tools / Articles pour travailler le plâtre


PLASTRIM

10' - 20.000 rpm/tpm

ø 1/10 mm
Pu / UE

		
	extracoarse, black ring/extragros, bague noire	
	070	070
	1	1
	61126	61124
Order-Nr.	●	●
Price in CAD	23.99	23.99


Application/Utilisation



MANDREL / MANDRIN

3' - 6.000 rpm/tpm

ø 1/10 mm
Pu / UE

	
	509077
	070
	12
Order-Nr.	●
Price in CAD	22.99

Follow us:



@dentologybyabutmentdirect

DENT**LOGY**
by **ABUTMENT DIRECT**
Same DNA. More Solutions.

**DENTOLOGY by
Abutment Direct Inc.**

16-250 Shields Court,
Markham, ON Canada L3R 9W7

Toll Free: 1-855-604-0465

Fax: (905) 604-0900

info@abutmentdirect.com

www.abutmentdirect.com/dentology